E.CROSS-A

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FR

NI.D1030A1-2025/05/06

EN

Instructions for use

IT

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ES

Nota de instrucciones

DE

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NL

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0. REVISION

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1. DESCRIPTION

> The E.CROSS -A cage is an intervertebral body fusion device intended to restore the intervertebral space and the physiological lordosis in skeletally mature patients. It is designed to replace lumbar discs and enable the fusion of adjacent vertebral body at vertebral levels L1 to S1. The cage must be packed with autograft or allograft material and must be implanted by healthcare professional.

> The E.CROSS-A intervertebral cage is intended to be inserted through an anterior or anterolateral (ALIF) approach. It must be used in combination with fixation screws or with a posterior or anterior fixation system.

> The cage is opened on both sides in order to contain the graft material and keep it in contact with the bone. The faces in contact with the vertebral endplates have a toothed shape to prevent any migration of the implant.

It is available in several heights, several lengths and several different shapes for various lordosis requirements.

> The E.CROSS-A cage is designed to be used with titanium bone fixing screws inserted into the vertebral bodies.

> The screws are available in several lengths.

> The cage is equipped with markers to facilitate the location and positioning of the implant using X-ray or other means of imaging.

> The materials used to manufacture this product satisfy current standard requirements (see the end of this leaflet).

2. INDICATIONS

The E.CROSS-A lumbar cage is indicated in the treatment of the following conditions of the lumbar spine (L1 to S1) in skeletally mature patients:

- Degenerative disc diseases and spinal instability at 1 or 2 levels;
- Spondylolisthesis (degenerative or isthmic lysis)

and retrolisthesis;

- Spinal stenosis;
- Revision;
- Pseudarthrosis.

3. CONTRAINDICATIONS

The list of contraindications includes:

- Severe osteoporosis;
- Skeletally immature;
- Prior fusion of treated segments;
- Acute or local infections, or allergies to any components;
- Unstable fractures (burst fractures, compression fractures), bone destruction;
- Spine tumor;
- Impossibility of reducing Spondylolisthesis to grade 1;
- Strong rotational deformation of pathological level.

4. WARNINGS / CAUTIONS

> The use of this product is reserved for surgeons with sufficient skills and experience in spinal surgery. The product must be used under the conditions specified by the manufacturer for the designated purpose.

It is highly recommended that the device be inserted by surgeons with substantial experience of the general problems associated with spine surgery and who have mastered surgical techniques using the anterior or lateral approach. It is incumbent on the user to adequately implement the proper resources to prevent neurological or vascular lesions that could result from the approach.

> Prior to use, surgeons shall be trained to the assembly and installation techniques and shall have read EUROS documentation. Insertion must be carried out in accordance with the recognised surgical technique. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly.

> Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and the selection and proper positioning of devices, are important factors in the successful use of this device or system by the surgeon. Furthermore, the selection of a suitable patient together with the

latter's cooperation are essential factors for a successful outcome.

Patients who present contra-indications must be turned down.

- > The surgeon should take into consideration factors that may affect the performance of the device or system:
 - The patient's weight, activity level and occupation;
 - The patient's bone state (osteoporosis, tumors);
 - The patient's diseases and disorders (infectious diseases, mental or neuromuscular disorders)
 - The patient's addiction (drug, alcoholism or nicotine);
 - The patient's pregnancy.

The surgeon must assess the patient's physical and mental state and decide whether the latter is suited to having the device fitted, can stand up to the operation and follow a suitable rehabilitation program.

- > The manufacturer is not liable for any complications arising from incorrect diagnosis, incorrect choice of device, the incorrect assembly of components of the device and/or incorrect operating techniques, the limitations of treatment methods or inadequate aseptic conditions.

- > Reporting an incident:

Any healthcare professional wishing to make a complaint or express their dissatisfaction in terms of the quality, identification, reliability, safety, efficacy and/or the performance of EUROS device must inform EUROS.

In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email (complaints@euros.fr). The competent authority of the Member State in which the user/or patient is established must be notified in case of a serious incident in relation to the device.

For any additional information involving the devices, please contact EUROS or your distributor.

EUROS SAS

Z.E. Athélia III - 13600 La Ciotat - France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

www.euros.fr

5. PRECAUTIONS

- > Caution should be exercised when handling the devices. The devices must not be scratched or damaged in any way whatsoever. Damaged devices and those with damaged packaging must not be used and should be returned to EUROS.
- > Devices must be kept in its original sealed packaging in a clean, dry location. Store the sterile products between 10 and 30°C. During transport, temperatures of up to 40°C for short periods can be tolerated. The packaging must not be directly exposed to the sun's rays or extreme temperatures.
- > Devices are for single use only and must not be reused. The reuse of a single use device that has come into contact with blood, bone, tissue or other body fluids may cause injury to the patient or user. The risks associated with the reuse of a single-use device may include, but are not limited to, mechanical failure and the transmission of infectious agents.
- > Devices must only be used with the instruments specified by EUROS in the reference surgical technique.
- > The E.CROSS-A cage is equipped with a space designed to accommodate an autogenous graft. The use of bone substitutes is subject to the responsibility of the surgeon.
- > Intervertebral fusion is effective over a maximum period of 2 years. The performance of E.CROSS-A cage is assured during this period of fusion. The mechanical strength of the cage has been checked for a two-year period after fitting. During this period, the grafted bone will have consolidated and taken over the mechanical function of the cage. However, if the cage needs to be removed after this period, every precaution should be taken to avoid damaging the surrounding structures during removal.
- > In the event of any disposal due to an error of storage or use, or to excision, devices should be disposed of following the disposal method for hospital waste in accordance with the procedures in force in the establishment.
- > Using all means available such as X-rays and templates, the surgeon must ensure that, by guaranteeing the concordance of scales between the X-rays and the templates, the implantation of the

device will cause no risks of impingement with other implants already in the patient's body.

6. EXPECTED PERFORMANCE AND BENEFITS

> Technical performances

- Enabling Implant position control (intraoperatively and postoperatively);
- Being easily inserted by anterior or anterolateral approach:
 - reliable connection with the implant holder;
 - Easy insertion design (rounded shape).
- Resisting to the mechanical loads applied over its lifetime:
 - Compression resistance ;
 - Shear resistance ;
 - Torsion resistance ;
 - Impaction resistance .
- Answering the implantation constraints
 - Biocompatible;
 - Being easily cleaned and sterilized ;
 - Being easily implantable and removable.

> Clinical performances

- Restoring the intervertebral space:
 - Being adapted to the different anatomies in height;
 - Being adapted to the different anatomies in lordosis.
- Enabling the vertebral bodies fusion:
 - Stability;
 - Optimal contact bone/graft surface ;
 - Avoid subsidence phenomenon ;
 - Assessment of the fusion by X-rays or other medical imaging.

> Clinical benefits

- Quality of life improvement, improving function scores ;
- Decreased pain;
- Restoration of spinal stability and physiological position.

7. ADVERSE EFFECTS / RESIDUAL RISKS

The list of potential adverse events includes :

- Migration of the implant
- Breakage of the implant

- Subsidence of the implant
- Patient's reaction to a foreign body
- Instability and degeneration of adjacent segments
- Pseudarthrosis
- Vertebral fracture
- Superficial and/or deep infection
- Injury to surrounding soft tissues (vascular, nervous, visceral, muscular)
- Dural leakage
- Neurological injury
- Paralysis
- Osteolysis
- Surgical complications (thrombosis, pulmonary embolism, myocardial infarction, bursitis, seroma, haematoma, difficulty healing)
- Pain
- Death
- *Intoxication, cellular damage, convulsion;
- *Cancer, tumor, immunological, and neurological disorders, genetic defects, non-heritable adverse effects in the progeny and/or an impairment of male or female reproductive functions or capacity.

Treatment of certain undesirable side effects could necessitate additional surgery.

**Note: these residual risks are related to the device biocompatibility. The biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable for the EUROS products.*

8. INFORMATION TO BE PROVIDED TO THE PATIENT

> The surgeon must inform the patient that the device's safety and durability will be affected by her/his weight, behaviour, namely the physical activities. Before leaving the care facility, the patient must receive instructions and warnings from the surgeon regarding the treatment and therapy to follow, and any restriction on physical and sporting activity or practices.

The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, as well as the need for regular medical follow-up.

> Implant card to be provide to the patient: Fill the

implant card following the associated instruction and give it to the patient. Explain him which information are available on the implant card. The surgeon must indicate to the patient that the additional information on the implant is available on the website by means of the product reference, in particular as regards to the precautions for use, the materials or MRI compatibility.

- > MR environment: the surgeon must inform the patient to indicate that he is wearing an implant before an MRI examination.

An MRI safety evaluation based on a bibliographical analysis and a demonstration of equivalence has demonstrated the E.CROSS-A cage is MR-conditional. A patient with this device can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla;
- Maximum spatial gradient field of 720-Gauss/cm;
- Total body average Specific Absorption Rate (SAR) maximum of 2 W/kg for 15 minutes of scanning (per pulse sequence).

MR imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the EUROS spine implant.

- > The surgeon must inform the patient to report any incident in relation to the devices to the surgeon.
- > The surgeon should inform the patient about the residual risks available in section 7 of these instructions for use.

9. CHOICE OF DEVICE / DEVICE COMPATIBILITY

- > Choosing the right cage size for each patient is a key element of this surgery. Once implanted, the cages are subjected to repeated efforts and their resistance is limited to their fit to the morphology of each patient.
- > The surgeon is responsible for the choice of implants, depending on the patient and the size of their vertebrae and on the correction sought by the surgeon. This choice may require preoperative planning. The instrumentation provided with the implant can also guide this choice thanks to a range of trial cages corresponding to the range of implants.
- > Information on device size is provided on the device

packaging and on the device itself. The surgical team must check this information before the surgery.

COMBINATION OF DEVICES

- > The surgeon must make sure that the inserted implant is compatible with devices already inserted in the patient.
- > The E.CROSS-A cage can be kept in place between the vertebrae without associated screws in the vertebral bodies by using an additional posterior or anterior stabilizing system.
- > It is recommended that the length of E.CROSS-A screw indicated on the cage's packaging label be used.

10. STERILITY

- > The E.CROSS-A cage is supplied sterile (gamma irradiation) at a sterility assurance level (SAL) of 10⁻⁶.
- > Inspect every package before use. Do not use the component if the integrity of the double envelope ensuring the maintenance of sterility has been compromised (loss of vacuum, damaged envelopes, damaged seal area), or if the expiry date has elapsed.
- > Once the package has been opened, the device must never be re-sterilised. Any re-sterilisation may have adverse effects on the device materials and cause premature defects.
- > For the additional fixings system and for instruments used to insert implants, see the relevant leaflet for information on sterilisation.

11. ADDITIONAL INFORMATION

The Summary of Safety and Clinical Performance is available upon request (email: qualite@euros.fr; Tel : +33 4 42 71 42 71).

Basic UDI-DI: 37004372R005BM

12. MEANING OF THE SYMBOLS

List of symbols applicable for the implants covered by these instructions for use:

	Sterilized using irradiation Double sterile barrier system
	Do not reuse

	Caution
	Batch code
	Keep dry
	Keep away from sunlight
	Storage temperature limitation
	Do not re-sterilise
	Use-by date
	Do not use if package is damaged and consult the instructions for use
	Catalogue number
	Manufacturer
	Date of manufacture
	Consult the instructions for use
	Medical device
	UDI
	Patient identification
	Name and address of the implanting healthcare institution and provider
	Date of implantation
	Information website for patients

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