# CROWN

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FR

NI.D1031A1-2024/10/04

EN

Instructions for use

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## INSTRUCTIONS FOR USE FOR SURGEONS AND MEDICAL STAFF

### 0. REVISION

Version & date	Modifications	URL internet
A1 - 04/10/2024	Creation	-

### 1. DESCRIPTION

- > The CROWN cranial vault is a custom-made device designed and manufactured from the patient's medical anatomical data for which it is intended. Under no circumstances may it be used for a recipient other than the patient for whom it was designed.
- > The CROWN device is intended to reconstruct the missing part of the cranial vault.
- > The CROWN cranial vault can be used for adult, adolescent or paediatric patients (with a fully developed cranial cavity).
- > The CROWN cranial vault has a constant thickness. It is bevelled on its contour to allow impingement-free positioning with the bone. A hole mesh is present on the entire surface to allow the surgeon to suspend the dura mater and reduce the vacuum spaces between the implant and the corpus callosum.
- > The materials used to manufacture these products satisfy current standard requirements (see the end of this leaflet).

### 2. INDICATIONS

The CROWN cranial vault is indicated for the reconstruction of cranial bone defects in the following cases:

- Trauma ;
- Tumour resection ;
- Decompressive craniectomy ;
- Revisions.

### 3. CONTRAINDICATIONS

The list of contraindications includes:

- Acute or chronic, local or systemic infections ;
- Missing cover tissues ;
- Allergy to any of the ingredients ;

- Unstable fractures, destruction, or poor bone quality that may affect the stability of the devices.

### 4. WARNINGS

- > The use of this product is reserved for surgeons with sufficient skills and experience in the cranioplasty field. The product must be used under the conditions specified by the manufacturer for the designated purpose.
- > Prior to use, surgeons shall be trained to the assembly techniques and shall have read EUROS documentation. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly.
- > The surgeon is responsible for ensuring that the imaging data provided as input data for the product design is consistent with the patient's anatomy at the time of implanting the custom device. The product design must be validated by surgeon before manufacturing. Imaging data of the patient are performed by trained operator following customization instruction.
- > Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and the selection and proper positioning of devices, are important factors in the successful use of this device or system by the surgeon. Furthermore, the selection of a suitable patient together with the latter's cooperation are essential factors for a successful outcome. Patients who present contra-indications must be turned down.
- > The surgeon should take into consideration factors that may affect the performance of the device:
  - The patient's activity level and occupation ;
  - The patient's bone state (osteoporosis, tumors) ;
  - The patient's allergies ;
  - The patient's diseases and disorders (infectious diseases, mental or neuromuscular disorders) ;
  - The patient's addiction (drug, alcoholism or nicotine).

The surgeon must assess the patient's physical and mental state and decide whether the latter is suited to having the device fitted, can stand up to the operation and follow a suitable rehabilitation program.

- > The manufacturer is not liable for any complications arising from incorrect diagnosis and/or incorrect operating techniques, the limitations of treatment methods or inadequate aseptic conditions.
- > MRI information: The safety and compatibility of the CROWN device have not been evaluated in a magnetic resonance (MR) environment. Heating, migration, or image artefacts for the CROWN device was not evaluated in a magnetic resonance (MR) environment.  
Nevertheless, the CROWN implant is manufactured from non-ferromagnetic materials (see the end of this leaflet).

#### > Reporting an incident:

Any healthcare professional wishing to make a complaint or express their dissatisfaction in terms of the quality, identification, reliability, safety, efficacy and/or the performance of EUROS device must inform EUROS.

In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email ([complaints@euros.fr](mailto:complaints@euros.fr)). The competent authority of the Member State in which the user/or patient is established must be notified in case of a serious incident in relation to the device.

For any additional information involving the devices, please contact EUROS or your distributor.

EUROS SAS

Z.E. Athélia III - 13600 La Ciotat - France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

[www.euros.fr](http://www.euros.fr)

## 5. PRECAUTIONS

- > Caution should be exercised when handling the devices. The devices must not be scratched or damaged in any way whatsoever. Damaged devices and those with damaged packaging must not be used and should be returned to EUROS.
- > Devices must be kept in its original sealed packaging in a clean, dry location. Store the sterile products between 10 and 30°C. During transport, temperatures of up to 40°C for short periods can be tolerated. The packaging must not be directly exposed to the sun's rays or extreme temperatures.

- > Devices are for single use only and must not be reused. The reuse of a single use device that has come into contact with blood, bone, tissue or other body fluids may cause injury to the patient or user. The risks associated with the reuse of a single-use device may include, but are not limited to, mechanical failure and the transmission of infectious agents.
- > In the event of any disposal due to an error of storage or use, or to excision, devices should be disposed of following the disposal method for hospital waste in accordance with the procedures in force in the establishment.
- > Using all means available such as X-rays and templates, the surgeon must ensure that, by guaranteeing the concordance of scales between the X-rays and the templates, the implantation of the device will cause no risks of impingement with other implants already in the patient's body. For that, the surgeon must remove the device previously implanted before implanting the new one.
- > A marking on the outer surface of the device indicates the orientation for implantation.

## 6. EXPECTED PERFORMANCE AND BENEFITS

### > Technical performances

- Allows reconstruction of a cranial vault bone defect:
  - Fills inside of the bone defect ;
  - Rebuilds the aesthetics of the outer table ;
  - Ensures the stability of the implant on the vault.
- Answering the implantation constraints
  - Biocompatible ;
  - Being easily cleaned and sterilized ;
  - Being easily implantable and removable.

### > Clinical performances

- Completes haemostasis and prevents extradural postoperative detachments :
  - Enables suspension of the dura mater.

### > Clinical benefits

- Reduced pain ;
- Aesthetics ;
- Improved quality of life.

## 7. ADVERSE EFFECTS / RESIDUAL RISKS

The list of potential adverse effects includes :

- Infection ;
- Migration of the implant ;
- Immune / allergic reactions ;
- Venous thrombosis, pulmonary embolism ;
- Cardiovascular disorders ;
- Haematoma ;
- Superficial wound / delayed wound healing ;
- Pain ;
- Discomfort ;
- Inflammatory reactions ;
- Neurological damage ;
- Soft tissues damage including vascular damage ;
- Autoimmune disease ;
- Formation of erythema or edema ;
- Implant rejection ;
- \*Intoxication, cellular damage, convulsion ;
- \*Cancer, tumor, immunological, and neurological disorders, genetic defects, non-heritable adverse effects in the progeny and/or an impairment of male or female reproductive functions or capacity.

Treatment of certain undesirable side effects could necessitate additional surgery.

*\*Note: these residual risks are related to the device biocompatibility. The biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable for the EUROS products.*

## 8. INFORMATION TO BE PROVIDED TO THE PATIENT

> The surgeon should inform the patient that the device's safety and durability will be affected by her/his behaviour, namely the physical activities. Before leaving the care facility, the patient must receive instructions and warnings from the surgeon regarding the treatment and therapy to follow, and any restriction on physical and sporting activity or practices. The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, but also of the need for regular medical follow-up. In particular, it is recommended not to perform physical activities which would lead to deterioration of the cranial vault and to use protective items accordingly. The surgeon should

inform the patient that it is not recommended to scratch the wound or scar after the procedure.

- > The lifespan of the implant may be shorter than the patient's lifetime. The surgeon must inform the patient that the device cannot and does not reproduce the flexibility, strength, reliability or durability of normal healthy bone, and that the device may need to be replaced in the future.
- > Implant card to be provide to the patient: fill the implant card using an indelible ink pen following the mentioned instruction and give it to the patient. Explain him which information are available on the implant card. The surgeon must indicate to the patient that the additional information on the implant is available on the website (accessible with the link shown on the implant card) by means of the product reference, in particular as regards to the precautions for use, the materials or MRI compatibility. The surgeon must inform the patient that they should always have their implant card with them.
- > The surgeon must inform the patient to report any incident in relation to the devices to the surgeon.
- > The surgeon should inform the patient about the residual risks available in section 7 of this instruction for use.

## 9. DEVICE COMPATIBILITY

- > The fixation of the cranial vault is the surgeon's choice. Most often, the fixation is made by screws and the plate used for cranio-maxillofacial surgery. Suspension of the dura mater can be performed with surgical wire.

## 10. STERILITY

- > The custom CROWN cranial vault is supplied sterile (gamma irradiation) at a sterility assurance level (SAL) of  $10^{-6}$ .
- > Inspect every package before use. Do not use the component if the integrity of the double envelope ensuring the maintenance of sterility has been compromised (loss of vacuum, damaged envelopes, damaged seal area), or if the expiry date has elapsed.
- > Once the package has been opened, the device must never be re-sterilised. Any re-sterilisation

may have adverse effects on the device materials and cause premature defects.


	Information website for patients
--	----------------------------------

## 11. ADDITIONAL INFORMATION

CROWN basic UDI-DI : 37004372C00184

## 12. MEANING OF THE SYMBOLS

List of symbols applicable for the devices covered by this instruction for use

	Sterilized using irradiation Double sterile barrier system
	Do not reuse
	Caution
	Batch code
	Keep dry
	Keep away from sunlight
	Storage temperature limitation
	Do not re-sterilise
	Use-by date
	Do not use if package is damaged and consult the instructions for use
	Catalogue number
	Manufacturer
	Date of manufacture
	Consult the instructions for use
	Medical device
	UDI
	Patient identification
	Name and address of the implanting healthcare institution and provider
	Date of implantation

## 13. MEANING OF THE ABBREVIATIONS

Qty: quantity - fr: French - en: English

## 14. MATERIAL

CROWN cranial vault	100% PEEK OPTIMA® LT1 (ASTM F2026)
------------------------	------------------------------------

**EUROS SAS**

Z.E. Athélia III - 13600 La Ciotat - France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

[www.euros.fr](http://www.euros.fr)