
EUROSTEM

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NI.D2005A1 – 2024/02/27

FR

EN

Instructions for use

IT

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ES

Nota de instrucciones

DE

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NL

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PT

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INSTRUCTIONS FOR USE
FOR SURGEONS AND MEDICAL STAFF

0. REVISION

| Version & date | Modifications | URL internet |
|--------------------|---------------|---------------------------------|
| A1 – 27/02/2024 | Creation | eifu.euros.fr/NI.D200 5A.pdf |

1. DESCRIPTION

> EUROS hip prostheses are intended to reconstruct the coxofemoral joint.

The femoral stems EUROSTEM are used for partial or total hip replacement.

The EUROSTEM femoral stem is available in five versions:

- Standard neck, cementless, collared ;
- Standard neck, cementless, collarless ;
- High-offset neck, cementless, collarless ;
- Coxa vara neck, cementless, collared ;
- Standard neck, cemented, collarless.

> The materials used to manufacture these products satisfy current standard requirements (see the end of this leaflet).

2. INDICATIONS

The indications for use of total hip replacement prostheses include:

- Severely painful and/or disabled joint as a result of osteoarthritis, post-traumatic arthritis, rheumatoid arthritis or congenital hip dysplasia ;
- Avascular necrosis of the femoral head ;
- Acute traumatic fracture of the femoral head or neck ;
- Failure of a previous hip surgery.

The indications for use of hemi-prostheses include:

- Acute traumatic fracture of the femoral head or neck ;
- Failure of a previous hip surgery.

3. CONTRAINDICATIONS

Hip replacement is contraindicated in the following cases:

- Acute, systemic or chronic infection ;
- Muscular, neurological or vascular deficiency of the affected limb ;
- Pathologies that may compromise the functionality of the implant in any way ;
- Bone destruction, or loss of bone characteristics that may compromise the stability of the implant ;

- Allergy to any prosthesis component ;
- For hemi-hip arthroplasty, any pathological condition of the acetabulum, that would preclude the use of the natural acetabulum as an appropriate articular surface for the hemi-hip prosthesis.

4. WARNINGS

> The use of this product is reserved for surgeons with sufficient skills and experience in hip replacement surgery. The product must be used under the conditions specified by the manufacturer for the designated purpose.

> Prior to use, surgeons shall be trained to the assembly and installation techniques and shall have read EUROS documentation. Insertion must be carried out in accordance with the recognised surgical technique. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly.

> Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and the selection and proper positioning of devices, are important factors in the successful use of this device or system by the surgeon. Furthermore, the selection of a suitable patient together with the latter's cooperation are essential factors for a successful outcome.

Patients who present contra-indications must be turned down.

> The surgeon should take into consideration factors that may affect the performance of the device or system:

- The patient's weight, activity level and occupation ;
- The patient's bone state (osteoporosis, tumors) ;
- The patient's allergies ;
- The patient's diseases and disorders (infectious diseases, mental or neuromuscular disorders) ;
- The patient's addiction (drug, alcoholism or nicotine).

The surgeon must assess the patient's physical and mental state and decide whether the latter is suited to having the device fitted, can stand up to the operation and follow a suitable rehabilitation program.

> The manufacturer is not liable for any complications arising from incorrect diagnosis, incorrect choice of device, the incorrect assembly of components of the device and/or incorrect operating techniques, the limitations of treatment methods or inadequate aseptic conditions.

> Reporting an incident:

Any healthcare professional wishing to make a complaint or express their dissatisfaction in terms of the

quality, identification, reliability, safety, efficacy and/or the performance of EUROS device must inform EUROS.

In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email (complaints@euros.fr). The competent authority of the Member State in which the user/or patient is established must be notified in case of a serious incident in relation to the device.

For any additional information involving the devices, please contact EUROS or your distributor.

EUROS SAS

Z.E. Athélia III - 13600 La Ciotat - France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

www.euros-orthopaedics.com

5. PRECAUTIONS

- > Caution should be exercised when handling the devices. The devices must not be scratched or damaged in any way whatsoever. Damaged devices and those with damaged packaging must not be used and should be returned to EUROS.
- > Devices must be kept in its original sealed packaging in a clean, dry location. Store the sterile products between 10 and 30°C. During transport, temperatures of up to 40°C for short periods can be tolerated. The packaging must not be directly exposed to the sun's rays or extreme temperatures.
- > Devices are for single use only and must not be reused. The reuse of a single use device that has come into contact with blood, bone, tissue or other body fluids may cause injury to the patient or user. The risks associated with the reuse of a single-use device may include, but are not limited to, mechanical failure and the transmission of infectious agents.
- > Devices must only be used with the instruments specified by EUROS in the reference surgical technique.
- > In the event of any disposal due to an error of storage or use, or to excision, devices should be disposed of following the disposal method for hospital waste in accordance with the procedures in force in the establishment.
- > Using all means available such as X-rays and templates, the surgeon must ensure that, by guaranteeing the concordance of scales between the X-rays and the templates, the implantation of the device will cause no risks of impingement with other implants already in the patient's body.

FEMORAL STEM

- > The EUROSTEM stems (size 8 and 9) are not recommended for patients over 55 kg.
- > Due to the metaphyseal anchorage of the stems, it is recommended to prepare the implant shoulder to prevent a wrong positioning.

FIXING THE BIOLOX DELTA HEAD TO THE FEMORAL STEM

- > The BIOLOX® delta head must be fixed very carefully to the EUROS femoral stem to ensure that the implant operates safely.
- > Test heads with the relevant diameter are used to check joint mobility and stability. Ensure that these test heads are removed.
- > Before fixing the BIOLOX® delta head to the EUROS femoral stem:
 - Thoroughly rinse the stem cone ;
 - Thoroughly dry the stem cone ;
 - Thoroughly inspect the stem cone and the female cone on the BIOLOX® delta head. Remove any foreign bodies such as residual tissue or bone or cement particles ;
 - The BIOLOX® delta head is placed manually on the stem cone, turning it slightly whilst pressing on the axis until it is firmly in position ;
 - Now place the plastic impactor on the pole of the BIOLOX® delta head and permanently attach the latter to the stem cone with a light stroke in the axial direction. The surface structure of the metallic cone undergoes plastic deformation on impact, which ensures optimal pressure distribution and good resistance to torsion on fixation.

ATTENTION :

Never tap the BIOLOX® delta head with a metallic impactor. Only use the designated plastic impactor supplied by EUROS.

6. EXPECTED PERFORMANCE AND BENEFITS

- > Technical performances
 - Enabling kinematics and hip joint mobility :
 - Restoring the femoral articulating surface ;
 - Being adapted to the different anatomies thanks to a large size range ;
 - Enabling joint stability :
 - Primary stability: stability of press-fit ;
 - Secondary stability: solidity and durability of cement or coating ;
 - Enabling reliable assembly with EUROS femoral

head (taper) :

- Matching perfectly with the cone of EUROS head prosthesis ;
- Torsion resistant ;
- Resisting to mechanical loads under normal conditions of use :
 - Endurance ;
- Answering the implantation constraints :
 - Biocompatible ;
 - Being easily cleaned and sterilized ;
 - Being easily implantable and removable.

> Clinical performances

- Restoration of hip articulation kinematics and mobility ;
- Restoration of hip joint stability.

> Clinical benefits

- Decreased pain ;
- Increased quality-of-life ;
- Easier walking.

7. ADVERSE EFFECTS / RESIDUAL RISKS

The list of potential adverse events associated with any hip arthroplasty includes :

- Excessive wear of the implant components secondary to impingement of components or damage of articular surfaces ;
- Instability ;
- Fracture, migration, loosening, dissociation, subsidence, protrusion, impingement, subluxation, or dislocation of the prosthesis or any of its components ;
- Possible detachment of the coating(s) on the femoral stem or acetabular shell components, potentially leading to increased debris particles ;
- Periprosthetic bone fracture, including femoral or acetabular perforation while seating the device ;
- Acetabular erosion or protrusion ;
- Osteolysis, geode ;
- Soft or surrounding tissues injuries ;
- Metal sensitivity reactions or other allergic/histological reactions to implant materials ;
- Inflammation, formation of erythema or oedema or implant rejection ;
- Superficial or deep infection, septic loosening ;
- Delayed wound healing ;
- Vascular damage resulting in blood loss and/or hematoma, potentially requiring transfusion ;

- Neurologic injury or neuropathy resulting in transient or permanent weakness, pain, and/or numbness ;
- Undesirable leg lengthening or shortening, claudication ;
- Periarticular calcification or heterotopic ossification, with or without impediment to joint mobility ;
- Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction ;
- Gastrointestinal or genitourinary complications ;
- Aggravation of other joint or back conditions due to positioning during surgery or neurological injury ;
- Decreased range of motion ;
- Squeaking (for ceramic-on-ceramic bearings) ;
- Intractable pain ;
- *Intoxication, cellular damage, convulsion ;
- *Cancer, tumor, immunological, and neurological disorders, genetic defects, non-heritable adverse effects in the progeny and/or an impairment of male or female reproductive functions or capacity.

Treatment of certain undesirable side effects could necessitate additional surgery.

**Note: these residual risks are related to the device biocompatibility. The biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable for the EUROS products.*

8. INFORMATION TO BE PROVIDED TO THE PATIENT

- > The surgeon should inform the patient that the device's safety and durability will be affected by her/his weight, behaviour, namely the physical activities. Before leaving the care facility, the patient must receive instructions and warnings from the surgeon regarding the treatment and therapy to follow, and any restriction on physical and sporting activity or practices. The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, but also of the need for regular medical follow-up. The lifespan of the implant may be shorter than the patient's lifetime or any particular given length of time. The surgeon must inform the patient that the device cannot and does not reproduce the flexibility, strength, reliability or durability of normal healthy bone, and that the device may need to be replaced in the future.
- > Implant card to be provide to the patient: fill the implant card using an indelible ink pen following the

mentioned instruction and give it to the patient. Explain him which information are available on the implant card. The surgeon must indicate to the patient that the additional information on the implant is available on the website (accessible with the link shown on the implant card) by means of the product reference, in particular as regards to the precautions for use, the materials or MRI compatibility. The surgeon must inform the patient that they should always have their implant card with them.

- > MR environment: the surgeon must inform the patient to indicate that he is wearing an implant before an MRI examination.

An MRI safety evaluation based on a bibliographical analysis and a demonstration of equivalence has demonstrated the EUROS hip prostheses are MR-conditional. A patient with these devices can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0-Tesla (3.0T);
- Maximum spatial gradient field of 1500-Gauss/cm ;
- Normal operating mode only ;
- Whole-body-averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning for patient in a 3T MR system.

MR imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the EUROS hip prosthesis.

- > The surgeon must inform the patient to report any incident in relation to the devices to the surgeon.
- > The surgeon should inform the patient about the residual risks available in section 7 of this instruction for use.

9. CHOICE OF DEVICE / DEVICE COMPATIBILITY

- > The surgeon is responsible for selecting the devices according to the patient's morphology and the nature of the condition.

This choice may require preoperative planning. The instrumentation provided with the implant can also guide this choice thanks to a range of trial implants corresponding to the range of implants.

- > Information on device size is provided on the device packaging and on the device itself. The surgical team must check this information before the surgery.
- > The use of a cemented or cementless implant depends on the surgeon's choice.

List of cemented/cementless implants of products:

| | Cemented | Cementless |
|----------|----------|------------|
| EUROSTEM | ✓ | ✓ |

- > In the case of devices designed to be implanted with cement or bone substitute, no incompatibility of use has been detected. The selection of the cement is the discretion of the surgeon.

Regarding devices designed to be implanted without cement, the manufacturer cannot, under any circumstances, be held responsible for implanting the device with cement.

COMBINATION OF DEVICES

- > EUROS devices must be combined only with EUROS devices. All combination with devices of other manufacturer is forbidden.
- > The combination head/stem is ensured by a 12/14 cone.

10. STERILITY

- > EUROS hip prostheses are supplied sterile (gamma irradiation) at a sterility assurance level (SAL) of 10⁻⁶.
- > Inspect every package before use. Do not use the component if the integrity of the double envelope ensuring the maintenance of sterility has been compromised (loss of vacuum, damaged envelopes, damaged seal area), or if the expiry date has elapsed.
- > Once the package has been opened, the device must never be re-sterilised. Any re-sterilisation may have adverse effects on the device materials and cause premature defects.
- > For instruments used to insert implants, see the relevant leaflet for information on sterilisation.

11. ADDITIONAL INFORMATION

The Summary of Safety and Clinical Performance is available upon request (email: qualite@euros.fr ; Tel : +33 4 42 71 42 71).

EUROSTEM basic UDI-DI : 37004372H0089M

12. MEANING OF THE SYMBOLS

List of symbols applicable for the devices covered by this instruction for use

| | |
|--|---|
|  | Sterilized using irradiation Double sterile barrier system |
|  | Do not reuse |

| | |
|--|--|
| | Caution |
| | Batch code |
| | Keep dry |
| | Keep away from sunlight |
| | Storage temperature limitation |
| | Do not re-sterilise |
| | Use-by date |
| | Do not use if package is damaged and consult the instructions for use |
| | Catalogue number |
| | Manufacturer |
| | Date of manufacture |
| | Consult the instructions for use |
| | Medical device |
| | UDI |
| | Patient identification |
| | Name and address of the implanting healthcare institution and provider |
| | Date of implantation |
| | Information website for patients |

13. MEANING OF THE ABBREVIATIONS

Qty: quantity - fr: French - en: English - it: Italian - es: Spanish - de: German – nl : Dutch – pt : Portuguese
 SS: Stainless steel - HAp: Hydroxyapatite powder – Ti: Titanium

14. MATERIAL

| | |
|---|--|
| EUROSTEM - STD - Collarless - Cemented | 100% Stainless steel X4CrNiMnMo 21-9-4 (ISO 5832-9 and ASTM F1586) |
| EUROSTEM- STD - Collarless - Cementless | 97.08 – 98.20% Ti6Al4V ELI (ISO 5832-3 and ASTM F620) + 1.80 – 2.92% HAp (ISO 13779-6) |

| | |
|--|--|
| EUROSTEM - STD - Collared - Cementless | 97.01 – 98.20% Ti6Al4V ELI (ISO 5832-3 and ASTM F620) + 1.80 – 2.99% HAp (ISO 13779-6) |
| EUROSTEM - High Offset - Collarless - Cementless | 97.20 – 98.28% Ti6Al4V ELI (ISO 5832-3 and ASTM F620) + 1.72 – 2.80% HAp (ISO 13779-6) |
| EUROSTEM – Coxa Vara - Collared - Cementless | 97.17 – 98.25% Ti6Al4V ELI (ISO 5832-3 and ASTM F620) + 1.75 – 2.83% HAp (ISO 13779-6) |

MATERIAL CHEMICAL COMPOSITION

Stainless steel X4CrNiMnMo (ISO 5832-9 and ASTM F1586)

| %C | %Si | %Mn | %P | %S | %N | %Cr | %Mo | %Ni | %Nb | %Cu | %Fe |
|-------|-------|-----------|--------|-------|-------------|------------|--------|---------|-------------|-------|---------|
| ≤0.08 | ≤0.75 | 2 to 4.25 | ≤0.025 | ≤0.01 | 0.25 to 0.5 | 19.5 to 22 | 2 to 3 | 9 to 11 | 0.25 to 0.8 | ≤0.25 | Balance |

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