# EUROSCUP FIXE & MOBILE

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FR

NI.D2007A2-2025/01/21

EN

Instructions for use

IT

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ES

Nota de instrucciones

DE

Gebrauchsanweinsug

NL

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PT

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## 0. REVISION

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## 1. DESCRIPTION

> EUROS hip prostheses are intended to reconstruct the coxofemoral joint.

The acetabular components are used for total hip replacement and exist in different versions:

- EUROSCUP FIXE cementless press-fit cup,
- EUROSCUP MOBILE cemented and cementless (Standard, Reconstruction and Tripode versions) dual mobility cup.

> The materials used to manufacture these products satisfy current standard requirements (see the end of this leaflet).

## 2. INDICATIONS

The indications for use of total hip replacement prostheses include:

- Severely painful and/or disabled joint as a result of osteoarthritis, post-traumatic arthritis, rheumatoid arthritis or congenital hip dysplasia
- Avascular necrosis of the femoral head,
- Acute traumatic fracture of the femoral head or neck,
- Failure of a previous hip surgery.

NB: The dual mobility design results in higher intra-prosthetic stability to address the treatment of patients with a high risk of dislocation (especially for elderly patients) or patients with recurrent dislocation.

## 3. CONTRAINDICATIONS

Hip replacement is contraindicated in the following cases:

- Acute, systemic or chronic infection,
- Muscular, neurological or vascular deficiency of the affected limb,
- Pathologies that may compromise the functionality of the implant in any way

- Bone destruction, or loss of bone characteristics that may compromise the stability of the implant,
- Allergy to any prosthesis component.

## 4. WARNINGS / CAUTIONS

> The use of this product is reserved for surgeons with sufficient skills and experience in hip replacement surgery. The product must be used under the conditions specified by the manufacturer for the designated purpose.

> Prior to use, surgeons shall be trained to the assembly and installation techniques and shall have read EUROS documentation. Insertion must be carried out in accordance with the recognised surgical technique. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly.

> Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and the selection and proper positioning of devices, are important factors in the successful use of this device or system by the surgeon. Furthermore, the selection of a suitable patient together with the latter's cooperation are essential factors for a successful outcome.

Patients who present contra-indications must be turned down.

> The surgeon should take into consideration factors that may affect the performance of the device or system:

- The patient's weight, activity level and occupation
- The patient's bone state (osteoporosis, tumors)
- The patient's allergies,
- The patient's diseases and disorders (infectious diseases, mental or neuromuscular disorders)
- The patient's addiction (drug, alcoholism or nicotine)

The surgeon must assess the patient's physical and mental state and decide whether the latter is suited to having the device fitted, can stand up to

the operation and follow a suitable rehabilitation program.

incorrect diagnosis, incorrect choice of device, the incorrect assembly of components of the device and/or incorrect operating techniques, the limitations of treatment methods or inadequate aseptic conditions.

> In very rare cases, a ceramic implant component may rupture. This may be due to prosthesis overload, e.g. due to the incorrect attachment of the BIOLOX<sup>®</sup> delta head on the stem cone, imperfect or incomplete adjustment between the BIOLOX<sup>®</sup> delta head and the stem cone, improper fixation of the BIOLOX<sup>®</sup> delta insert in the metal cup or lack of attention to detail in terms of the recommended positioning of BIOLOX<sup>®</sup> delta inserts (inclination/anteversion).

The use of prosthetic components not supplied by EUROS may also cause the BIOLOX<sup>®</sup> delta head or insert to fracture. To minimise this risk, all BIOLOX<sup>®</sup> delta heads and inserts are individually inspected prior to delivery.

Substantial, temporary overloads in the event of falls, accidents or excessive overload due to an intense sporting activity, for instance, can cause the implant to rupture, often long after the event in question. If a ceramic component is broken, the matching of metal (head) and polyethylene (insert) as well as metal with metal is contraindicated during revision.

> Reporting an incident:

Any healthcare professional wishing to make a complaint or express their dissatisfaction in terms of the quality, identification, reliability, safety, efficacy and/or the performance of EUROS device must inform EUROS.

In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email ([complaints@euros.fr](mailto:complaints@euros.fr)). The competent authority of the Member State in which the user/or patient is established must be notified in case of a serious incident in relation to the device.

For any additional information involving the devices, please contact EUROS or your distributor.

> The manufacturer is not liable for any complications arising from

EUROS SAS

Z.E. Athélia III - 13600 La Ciotat - France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

[www.euros.fr](http://www.euros.fr)

## 5. PRECAUTIONS

> Caution should be exercised when handling the devices. The devices must not be scratched or damaged in any way whatsoever. Damaged devices and those with damaged packaging must not be used and should be returned to EUROS.

> The BIOLOX<sup>®</sup> delta heads and inserts are extremely sensitive to degradation.

Scratches or shock may cause excessive wear and tear or rupture, and therefore cause complications. Particularly cautious handling is therefore indicated.

> Devices must be kept in its original sealed packaging in a clean, dry location. Store the sterile products between 10 and 30°C. During transport, temperatures of up to 40°C for short periods can be tolerated. The packaging must not be directly exposed to the sun's rays or extreme temperatures.

> Devices are for single use only and must not be reused. The reuse of a single use device that has come into contact with blood, bone, tissue or other body fluids may cause injury to the patient or user. The risks associated with the reuse of a single-use device may include, but are not limited to, mechanical failure and the transmission of infectious agents.

> A BIOLOX<sup>®</sup> delta insert or head that has already been used must not be reused.

Previously used BIOLOX<sup>®</sup> delta heads and inserts may have defects that are not visible to the naked eye. Defects of any kind can have a negative impact on the functionality and/or stability of the ceramic component. Therefore, safety in use cannot be guaranteed. Defects (e.g. impact points or metal deposition) can cause excessive wear and tear or fracture and lead to complications.

Consequently, only BIOLOX® delta inserts and heads removed from their original packaging immediately prior to implantation can be used.

Given the precise necessary adjustment between the BIOLOX delta insert and the femoral stem, only new intact components can be combined. This also means, for example, that a BIOLOX® delta head that has been placed on a stem and then removed should not be placed on another stem. Similarly, a BIOLOX® delta insert that has been damaged must not be used and should be disposed of. This also applies to a BIOLOX® delta head or insert that has fallen on the floor, for instance.

- > Devices must only be used with the instruments specified by EUROS in the reference surgical technique.
- > In the event of any disposal due to an error of storage or use, or to excision, devices should be disposed of following the disposal method for hospital waste in accordance with the procedures in force in the establishment.
- > Using all means available such as X-rays and templates, the surgeon must ensure that, by guaranteeing the concordance of scales between the X-rays and the templates, the implantation of the device will cause no risks of impingement with other implants already in the patient's body.

### **CEMENTLESS CUPS**

- > The metal cup provides a press-fit versus the last acetabular reamer used. Hence, when dealing with sclerotic bone, all possible precautions should be taken to prevent bone fracture.

### **EUROSCUP FIXE**

- > Primary fixation can be reinforced by Ø6,5 optional cancellous bone screws.

### **EUROSCUP MOBILE**

- > Primary fixation can be reinforced by impacted pegs (EUROSCUP MOBILE Tripode) and/or Ø4,5 cancellous screws through iliac flange(s) (EUROSCUP MOBILE Tripode and Reconstruction versions).

### **FIXING BIOLOX® delta INSERTS IN THE EUROSCUP FIXE CUP**

- > The BIOLOX® delta insert must be fixed very carefully in the acetabular cup to ensure that the device operates safely.

- > Before making a final decision regarding the implantation of a ceramic insert, a plastic test insert will initially be placed in the acetabulum and the components of the implanted stem. Ensure that this insert is then removed.
- > Before placing the BIOLOX® delta insert in the EUROSCUP FIXE cup, thoroughly rinse the attachment area and dry carefully. Then remove any foreign bodies from the prosthetic components such as tissue residues, bone or cement particles.
- > Pay particular attention when centring the BIOLOX® delta insert in EUROSCUP FIXE cup. Follow the EUROSCUP FIXE operating instructions.

### **ATTENTION :**

Never strike a BIOLOX® delta insert with a metallic hammer. Only use the designated plastic impactor supplied by EUROS.

## **6. EXPECTED PERFORMANCE AND BENEFITS**

- > Technical performances
  - Enabling kinematics and hip joint mobility:
    - Restoring the acetabular articulating surface.
  - Enabling joint stability:
    - Primary stability: stability of press-fit;
    - Secondary stability: solidity and durability of cement or coating or additional screw fixation.
  - Enabling adequate joint tension:
    - Being adapted to the different anatomies with external diameter variation;
  - Enabling reliable assembly with EUROS femoral heads:
    - Matching perfectly with the femoral head shape;
    - Dislocation resistant.
  - Resisting to wear under normal conditions of use.
  - Answering the implantation constraints :
    - Biocompatible;
    - Being easily cleaned and sterilized;
    - Being easily implantable and removable;

- Being visible through medical imaging;
  - Being MRI compatible.
- > Clinical performances
- Restoration of hip articulation kinematics and mobility;
  - Restoration of hip joint stability.
- > Clinical benefits
- Decreased pain ;
  - Increased quality-of-life ;
  - Easier walking.

## 7. ADVERSE EFFECTS / RESIDUAL RISKS

The list of potential adverse events associated with any hip arthroplasty:

- Excessive wear of the implant components secondary to impingement of components or damage of articular surfaces;
- Instability;
- Fracture, migration, loosening, dissociation, subsidence, protrusion, impingement, subluxation, or dislocation of the prosthesis or any of its components;
- Possible detachment of the coating(s) on the femoral stem or acetabular shell components, potentially leading to increased debris particles;
- Periprosthetic bone fracture, including femoral or acetabular perforation while seating the device;
- Acetabular erosion or protrusion;
- Osteolysis, geode;
- Soft or surrounding tissues injuries;
- Metal sensitivity reactions or other allergic/histological reactions to implant materials; inflammation, formation of erythema or oedema or implant rejection;
- Superficial or deep infection, septic loosening;
- Delayed wound healing;
- Vascular damage resulting in blood loss and/or hematoma, potentially requiring transfusion;
- Neurologic injury or neuropathy resulting in transient or permanent weakness, pain, and/or numbness;
- Muscle or soft tissues injury;
- Undesirable leg lengthening or shortening,

- claudication;
- Periarticular calcification or heterotopic ossification, with or without impediment to joint mobility;
- Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction;
- Gastrointestinal or genitourinary complications;
- Aggravation of other joint or back conditions due to positioning during surgery or neurological injury;
- Decreased range of motion;
- Squeaking (for ceramic-on-ceramic bearings);
- Intractable pain;
- \*Intoxication, cellular damage, convulsion
- \*Cancer, tumor, immunological, and neurological disorders, genetic defects, non-heritable adverse effects in the progeny and/or an impairment of male or female reproductive functions or capacity

Treatment of certain undesirable side effects could necessitate additional surgery.

*\*Note: these residual risks are related to the device biocompatibility. The biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable for the EUROS products.*

## 8. INFORMATION TO BE PROVIDED TO THE PATIENT

- > The surgeon should inform the patient that the device's safety and durability will be affected by her/his weight, behaviour, namely the physical activities. Before leaving the care facility, the patient must receive instructions and warnings from the surgeon regarding the treatment and therapy to follow, and any restriction on physical and sporting activity or practices. The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, but also of the need for regular medical follow-up. The lifespan of the implant may be shorter than the patient's lifetime or any particular given length of time. The surgeon must

inform the patient that the device cannot and does not reproduce the flexibility, strength, reliability or durability of normal healthy bone, and that the device may need to be replaced in the future.

> Implant card to be provide to the patient: fill the implant card using an indelible ink pen following the mentioned instruction and give it to the patient. Explain him which information are available on the implant card. The surgeon must indicate to the patient that the additional information on the implant is available on the website (accessible with the link shown on the implant card) by means of the product reference, in particular as regards to the precautions for use, the materials or MRI compatibility. The surgeon must inform the patient that they should always have their implant card with them.

> MR environment: the surgeon must inform the patient to indicate that he is wearing an implant before an MRI examination.

An MRI safety evaluation based on a bibliographical analysis and a demonstration of equivalence has demonstrated the EUROS hip prostheses are MR-conditional. A patient with these devices can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0-Tesla
- Maximum spatial gradient field of 1500-Gauss/cm
- Normal operating mode only
- Whole-body-averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning for patient in a 3T MR system.

MR imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the EUROS hip prosthesis.

> The surgeon must inform the patient to report any incident in relation to the devices to the surgeon.

> CMR substances: the surgeon must inform the patient that the components made in CoCr alloy contain more than 0.1% (weight/weight) of cobalt identified as a carcinogenic and toxic for reproduction substance according to the regulation 1272/2008, 1907/2006 and 2020/2096. However, the cobalt chromium alloy has long clinical history and remains the best alternative in hip prosthesis. Our biological risk assessment

report demonstrated that the biological risk was controlled and considered acceptable.

> The surgeon should inform the patient about the residual risks available in section 7 of this instruction for use.

## 9. CHOICE OF DEVICE / DEVICE COMPATIBILITY

> The surgeon is responsible for selecting the devices according to the patient's morphology and the nature of the condition.

This choice may require preoperative planning. The instrumentation provided with the implant can also guide this choice thanks to a range of trial implants corresponding to the range of implants.

> Information on device size is provided on the device packaging and on the device itself. The surgical team must check this information before the surgery.

> The use of a cemented or cementless implant depends on the surgeon's choice.

> List of cemented/cementless implants of products:

		Cemented	Cementless
EUROSCUP FIXE	Cup		✓
EUROSCUP MOBILE	Cup	✓	✓
	EUROSCUP MOBILE tripode, impacted pegs		✓

> In the case of devices designed to be implanted with cement or bone substitute, no incompatibility of use has been detected. The selection of the cement is the discretion of the surgeon.

Regarding devices designed to be implanted without cement, the manufacturer cannot, under any circumstances, be held responsible for implanting the device with cement.

## COMBINATION OF DEVICES

> EUROS devices must only be combined with EUROS devices or with other devices distributed by EUROS.

- > Only inserts and cups of the same size and the same range can be combined.  
The BIOLOX® Delta inserts are compatible with the EUROSCUP FIXE cups.
- > The head / insert assembly must be made between head and insert of corresponding diameter.  
EUROSCUP FIXE combination between PE insert and femoral head (metallic or ceramic):

Femoral head diameter	PE insert size					
	42	44	46/48 /50	52/54	56/ 58	60/62/ 64/66
∅22 (metallic)	x	x				
∅28 (metallic or ceramic)			x	x	x	x
∅32 (metallic or ceramic)				x	x	x

EUROSCUP FIXE combination between ceramic insert and ceramic femoral head:

Femoral head diameter	Ceramic insert size			
	46/48/ 50	52/54	56/58	60/62/ 64/66
∅28				
∅32	x			
∅36		x	x	x
∅40			x	x

The BIOLOX® Delta inserts must only be used with BIOLOX® Delta femoral heads.

EUROSCUP MOBILE combination between PE insert and femoral head (metallic or ceramic):

Femoral head diameter	PE insert size												
	42	44	46	48	50	52	54	56	58	60	62	64	66
∅22 (metallic)	x	x	x	x	x								
∅28 (metallic or ceramic)			x	x	x	x	x	x	x	x	x	x	x

## 10. STERILITY

- > EUROS hip prostheses are supplied sterile (gamma irradiation) at a sterility assurance level (SAL) of 10<sup>-6</sup>.
- > Inspect every package before use. Do not use the component if the integrity of the double envelope ensuring the maintenance of sterility has been

compromised (loss of vacuum, damaged envelopes, damaged seal area), or if the expiry date has elapsed.

- > Once the package has been opened, the device must never be re-sterilised. Any re-sterilisation may have adverse effects on the device materials and cause premature defects.
- > For instruments used to insert implants, see the relevant leaflet for information on sterilisation.




## 11. ADDITIONAL INFORMATION

The Summary of Safety and Clinical Performance is available upon request (email: [qualite@euros.fr](mailto:qualite@euros.fr); Tel : +33 4 42 71 42 71).

EUROSCUP FIXE & MOBILE basic UDI-DI : 37004372H0049D

## 12. MEANING OF THE SYMBOLS

List of symbols applicable for the implants covered by these instructions for use

	Sterilized using irradiation Double sterile barrier system
	Do not reuse
	Caution
	Batch code
	Keep dry
	Keep away from sunlight
	Storage temperature limitation
	Do not re-sterilise
	Use-by date
	Do not use if package is damaged and consult the instructions for use
	Catalogue number
	Manufacturer

	Date of manufacture
	Consult the instructions for use
	Carcinogenic, mutagenic, toxic for reproduction or endocrine disruptor
	Medical device
	UDI
	Patient identification
	Name and address of the implanting healthcare institution and provider
	Date of implantation
	Information website for patients

### 13. MEANING OF THE ABBREVIATIONS

Qty: quantity - fr: French - en: English - it: Italian - es: Spanish - de: German - nl : Dutch - pt: Portuguese

CoCr: Cobalt chromium - HAp: Hydroxyapatite powder - Ti: Titanium - UHMWPE: ultra high molecular weight polyethylene - Ti6Al4V ELI : aluminum-vanadium alloy titanium

### 14. MATERIAL

EUROSCUP FIXE - Acetabular cup	86,48% - 96,09% Ti6Al4V ELI (ISO 5832-3 and ASTM F136) + 2,35% - 8,75% Ti coating (ISO 13179-1) + 1,5% - 5,15% HAp coating (ISO 13779-6)
Standard Insert	100% UHMWPE (ISO 5834-2 and ASTM F648)
Ceramic Insert	100% Ceramic (BILOX <sup>®</sup> delta)

EUROSCUP MOBILE - Acetabular cemented Cup	100% CoCr (ISO 5832-4 and ASTM F75)
EUROSCUP MOBILE - Acetabular cementless Cup	86,20% - 92,18% CoCr (ISO 5832-4 and ASTM F75) + 4,66% - 8,93% Ti coating (ISO 13179-1) + 2,94% - 5,26% HAp coating (ISO 13779-6)

EUROSCUP MOBILE - Reconstruction Acetabular Cup	87,6% - 92,86% CoCr (ISO 5832-4 and ASTM F75) + 4,27% - 8,01% Ti coating (ISO 13179-1) + 2,69% - 4,71% HAp coating (ISO 13779-6)
EUROSCUP MOBILE - Tripode Acetabular Cup	82,48% - 90,42% CoCr (ISO 5832-4 and ASTM F75) + 1,86% - 5,10% Ti6Al4V ELI (ISO 5832-3 and ASTM F136) + 4,51% - 8,28% Ti coating (ISO 13179-1) + 2,89% - 5,03% HAp coating (ISO 13779-6)
Dual Mobility Insert	100% UHMWPE (ISO 5834-2 and ASTM F648)

### MATERIAL CHEMICAL COMPOSITION

CoCr casting alloy (ISO 5832-4 and ASTM F75) [%m/m]

Cr	Mo	Ni	Fe	C	Si	Mn	W
27 to 30	5 to 7	≤0.5	≤0.75	≤0.35	≤1	≤1	≤0.2
P	S	N	Al	Ti	B	Co	P
≤0,020	≤0,010	≤0,25	≤0,10	≤0,10	≤0,010	Balance	≤0,020

UHMWPE (ISO 5834-2 and ASTM F648) [mg/kg]

Ash	Ti	Ca	Cl	Al
≤125	≤40	≤5	≤30	≤20

Ti6Al4V ELI (ISO 5832-3 and ASTM F136) [%m/m]

C	Al	V	H	O	N	Fe	Ti
≤0.08	5.5 to 6.5	3.5 to 4.5	≤0.012	≤0.13	≤0.05	≤0.25	Balance

HAp (ISO 13779-6) [mg/kg]

Pb	Hg	Bi	As	Sb	Sn
≤30	≤5	/	≤3	/	/
Cd	Ag	Cu	Mo	Total heavy metals	
≤5	/	/	/	≤30	

Ti coating (ISO 13179-1) [%m/m]

%C	%H	%Fe	%N	%O	%Ti
≤0.1	≤0.2	≤0.6	≤5	≤10	Balance

EC mark granted date: Ceramic insert (2014), EUROSCUP FIXE & MOBILE (2020).

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