EUROP UNI FIXE

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FR

NI.D3006A1 – 2023/12/15

EN

Instructions for use

IT

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ES

Nota de instrucciones

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INSTRUCTIONS FOR USE
FOR SURGEONS AND MEDICAL STAFF

0. REVISION

Version & date	Modifications	URL internet
A1 – 15/12/2023	Creation	eifu.euros.fr/NI.D300 6A.pdf

1. DESCRIPTION

- > The EUROP UNI FIXE knee joint prosthesis is intended to replace the articulating surface of the damaged/degenerative knee joint (unicompartmental knee replacement) in skeletally mature patients, to reduce pain, improve quality of life and knee mobility. It can be implanted with or without bone cement. The cemented version is intended to be used with bone cement. It shall be implanted by healthcare professional.
- > The materials used to manufacture these products satisfy current standard requirements (see the end of this leaflet).

2. INDICATIONS

EUROP UNI FIXE is indicated for primary partial (medial or lateral compartment) knee replacement surgery in skeletally mature patients with the following indications:

- Degenerative arthritis, or post-traumatic arthritis ;
- Moderate varus, valgus, or flexion deformity in which the ligamentous structures can be returned to adequate function and stability ;
- Avascular necrosis of the knee ;
- Post-traumatic loss of knee joint configuration and function.

In all cases, the joint cartilage of the opposite compartment must be intact and able to support normal loads.

3. CONTRAINDICATIONS

EUROP UNI FIXE is contraindicated in the following cases:

- Progressive, acute or chronic infections, local to the operative site or systemic, that may affect the prosthetic joint ;
- Any loss of musculature or serious lesions of muscles, nerves and/or blood vessels deficiency, which put the affected limb at risk ;
- Lacking bone substance or inadequate bone quality on femoral or tibial surfaces, which might compromise the stability of the prosthesis

component ;

- Skeletal immaturity ;
- Pathologies that may compromise the functionality of the implant in any way ;
- Inflammatory arthritis, rheumatoid arthritis ;
- Significant chondrocalcinosis of the articular surface ;
- Advanced involvement of the other compartments or a significant extra-articular deformity ;
- Over 30° of fixed varus or valgus deformity ;
- Ligamentous instability of the knee, failure of one or both cruciate ligaments ;
- Poor skin coverage around the knee that would make the procedure unjustifiable ;
- Neuropathic arthropathy (Charcot disease) ;
- Mental or neuromuscular disorders (which may create an unacceptable risk to the patient and can be a source of postoperative complications) ;
- Hypersensitivities or allergies to the materials used.

4. WARNINGS

- > The use of this product is reserved for surgeons with sufficient skills and experience in knee replacement surgery. The product must be used under the conditions specified by the manufacturer for the designated purpose.
- > Prior to use, surgeons shall be trained to the assembly and installation techniques and shall have read EUROS documentation. Insertion must be carried out in accordance with the recognised surgical technique. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly.
- > Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and the selection and proper positioning of devices, are important factors in the successful use of this device or system by the surgeon. Furthermore, the selection of a suitable patient together with the latter's cooperation are essential factors for a successful outcome.
- > Patients who present contra-indications must be turned down.
- > The surgeon should take into consideration factors that may affect the performance of the device or system:
 - The patient's weight, activity level and occupation ;
 - The patient's bone state (osteoporosis, tumors) ;
 - The patient's allergies ;
 - The patient's diseases and disorders (infectious diseases, mental or neuromuscular disorders) ;
 - The patient's addiction (drug, alcoholism or nicotine).

The surgeon must assess the patient's physical and mental state and decide whether the latter is suited to having the device fitted, can stand up to the operation and follow a suitable rehabilitation program.

- > The manufacturer is not liable for any complications arising from incorrect diagnosis, incorrect choice of device, the incorrect assembly of components of the device and/or incorrect operating techniques, the limitations of treatment methods or inadequate aseptic conditions.

> Reporting an incident:

Any healthcare professional wishing to make a complaint or express their dissatisfaction in terms of the quality, identification, reliability, safety, efficacy and/or the performance of EUROS device must inform EUROS. In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email (complaints@euros.fr). The competent authority of the Member State in which the user/or patient is established must be notified in case of a serious incident in relation to the device.

For any additional information involving the devices, please contact EUROS or your distributor.

EUROS SAS

Z.E. Athélia III - 13600 La Ciotat - France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

www.euros-orthopaedics.com

5. PRECAUTIONS

- > Caution should be exercised when handling the devices. The devices must not be scratched or damaged in any way whatsoever. Damaged devices and those with damaged packaging must not be used and should be returned to EUROS.
- > Devices must be kept in its original sealed packaging in a clean, dry location. Store the sterile products between 10 and 30°C. During transport, temperatures of up to 40°C for short periods can be tolerated. The packaging must not be directly exposed to the sun's rays or extreme temperatures.
- > Devices are for single use only and must not be reused. The reuse of a single use device that has come into contact with blood, bone, tissue or other body fluids may cause injury to the patient or user. The risks associated with the reuse of a single-use device may include, but are not limited to, mechanical failure and the transmission of infectious agents.
- > Devices must only be used with the instruments

specified by EUROS in the reference surgical technique.

- > In the event of any disposal due to an error of storage or use, or to excision, devices should be disposed of following the disposal method for hospital waste in accordance with the procedures in force in the establishment.
- > Using all means available such as X-rays and templates, the surgeon must ensure that, by guaranteeing the concordance of scales between the X-rays and the templates, the implantation of the device will cause no risks of impingement with other implants already in the patient's body.

Tibial baseplate

- > It is important to start with the tibial cut for two main reasons:
 - because the distal femoral cut is tibia dependent ;
 - to respect implants combination.
- > Avoid under sizing the tibial baseplate to avoid any risk of depression. The ideal solution is to aim for cortical support throughout the periphery.
- > Do not damage the cruciate ligaments when making the sagittal cut.
- > The tibial slope should be respected.
- > Securing with screws is strongly recommended when using an implant without cement.
- > Maintain strong pressure against the «Full PE» base plate until the cement has completely set.

Tibial insert

- > Position the tibial insert after placing the femoral condyle.

Femoral condyle

- > The distal femoral cut is tibia dependent.
- > Avoid hypercorrecting as this causes premature wear of the contralateral compartment.
- > Make sure that the top of the implant does not conflict with the patella.

6. EXPECTED PERFORMANCE AND BENEFITS

> Technical performances

- Restoring kinematics and knee joint mobility :
 - Respecting physiological ligament tension ;
 - Being adapted to the different anatomies thanks to a large size range ;
- Resistance to mechanical loads under normal conditions of use :
 - Compression resistance ;
 - Integrity of connections between components ;

- Resistance to wear under normal conditions of use ;
- Answering the implantation constraints :
 - Biocompatible ;
 - Being easily cleaned and sterilized ;
 - Being easily implantable and removable.

> Clinical performances :

- Improving knee mobility (depending on the indication) ;
- Reducing pain level during daily activities.

> Clinical benefits :

- Decreased pain ;
- Increased quality of life ;
- Easier walking.

7. ADVERSE EFFECTS / RESIDUAL RISKS

The list of potential adverse events associated with any knee arthroplasty includes :

- Loosening or breakage/damage of knee implant components ;
- Reaction of/damage to the tissues surrounding the prosthesis, blood loss ;
- Dislocation, luxation and/or instability of the joint ;
- Incorrect alignment or malposition of implant components (hypercorrection...) ;
- Bone fracture or nervous disorders ;
- Superficial wounds, hematoma, swelling, infection or inflammation ;
- Leg length discrepancy ;
- Low range of motion ;
- Pain ;
- Excessive implant wear or bone osteolysis ;
- Cardiovascular disorders, venous thromboembolic disease or pulmonary embolism ;
- Allergy or sensitivity to material ;
- Osteophytes formation ;
- Osteoarthritis progression ;
- *Intoxication, cellular damage, convulsion ;
- *Cancer, tumor, immunological, and neurological disorders, genetic defects, non-heritable adverse effects in the progeny and/or an impairment of male or female reproductive functions or capacity.

Treatment of certain undesirable side effects could necessitate additional surgery.

**Note: these residual risks are related to the device biocompatibility. The biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable for the EUROS products.*

8. INFORMATION TO BE PROVIDED TO THE PATIENT

> The surgeon should inform the patient that the device's safety and durability will be affected by her/his weight, behaviour, namely the physical activities. Before leaving the care facility, the patient must receive instructions and warnings from the surgeon regarding the treatment and therapy to follow, and any restriction on physical and sporting activity or practices. The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, but also of the need for regular medical follow-up. The lifespan of the implant may be shorter than the patient's lifetime or any particular given length of time. The surgeon must inform the patient and the parent(s) (or the parent(s)) that the device cannot and does not reproduce the flexibility, strength, reliability or durability of normal healthy bone, and that the device may need to be replaced in the future.

> Implant card to be provided to the patient: fill the implant card using an indelible ink pen following the mentioned instruction and give it to the patient. Explain him which information are available on the implant card. The surgeon must indicate to the patient that the additional information on the implant is available on the website (accessible with the link shown on the implant card) by means of the product reference, in particular as regards to the precautions for use, the materials or MRI compatibility. The surgeon must inform the patient that they should always have their implant card with them.

> MR environment: the surgeon must inform the patient to indicate that he is wearing an implant before an MRI examination.

An MRI safety evaluation based on a bibliographical analysis and a demonstration of equivalence has demonstrated the EUROS knee prostheses are MR-conditional. A patient with these devices can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 3.0-Tesla (3.0T) or less ;
- Maximum spatial gradient field of 750-Gauss/cm ;
- Whole-body-averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning for patient in a 3T MR system.

MR imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the EUROS knee device.

> The surgeon must inform the patient to report any incident in relation to the devices to the surgeon.

- > CMR substances: the surgeon must inform the patient that the baseplate and the condyle of EUROP UNI FIXE contain more than 0.1% (weight/weight) of cobalt identified as a carcinogenic and toxic for reproduction substance according to the regulation 1272/2008, 1907/2006 and 2020/2096. However, the cobalt chromium alloy has long clinical history and remains the best alternative in knee prosthesis. Our biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable.
- > The surgeon should inform the patient about the residual risks available in section 7 of this instruction for use.

9. CHOICE OF DEVICE / DEVICE COMPATIBILITY

- > The surgeon is responsible for selecting the devices according to the patient's morphology and the nature of the condition. This choice may require preoperative planning. The instrumentation provided with the implant can also guide this choice thanks to a range of trial implants corresponding to the range of implants.
- > Information on device size is provided on the device packaging and on the device itself. The surgical team must check this information before the surgery.
- > The use of a cemented or cementless implant depends on the surgeon's choice.

List of cemented/cementless implants of products:

	Cemented	Cementless
Condyle	✓	✓
Baseplate	✓	✓
Full PE Insert	✓	

- > In the case of devices designed to be implanted with cement or bone substitute, no incompatibility of use has been detected. The selection of the cement is the discretion of the surgeon.

Regarding devices designed to be implanted without cement, the manufacturer cannot, under any circumstances, be held responsible for implanting the device with cement.

- > We recommended you adjust the size of the tibial baseplate to avoid any risk of conflict with the lateral ligaments in the case of over-sizing, or subsidence of the baseplate in the case of under-sizing. The ideal size corresponds to cortical support on the entire periphery.

COMBINATION OF DEVICES

- > EUROS devices must be combined only with EUROS devices. All combination with devices of other manufacturer is forbidden.
- > Compatibility table representing the authorised combinations between the components for correct assembly:

		Femoral component size						
		1	2	3	4	5	6	7
Baseplate / Insert or Full PE	1	✓	✓					
	2	✓	✓	✓				
	3	✓	✓	✓	✓	✓	✓	✓
	4		✓	✓	✓	✓	✓	✓
	5			✓	✓	✓	✓	✓
	6					✓	✓	✓

10. STERILITY

- > EUROS devices are supplied sterile (gamma irradiation) at a sterility assurance level (SAL) of 10^{-6} .
- > Inspect every package before use. Do not use the component if the integrity of the double envelope ensuring the maintenance of sterility has been compromised (loss of vacuum, damaged envelopes, damaged seal area), or if the expiry date has elapsed.
- > Once the package has been opened, the device must never be re-sterilised. Any re-sterilisation may have adverse effects on the device materials and cause premature defects.
- > For instruments used to insert implants, see the relevant leaflet for information on sterilisation.

11. ADDITIONAL INFORMATION

The Summary of Safety and Clinical Performance is available upon request (email: qualite@euros.fr ; Tel : +33 4 42 71 42 71).

EUROP UNI FIXE basic UDI-DI : 37004372G00292

12. MEANING OF THE SYMBOLS

List of symbols applicable for the devices covered by this instruction for use

	Sterilized using irradiation Double sterile barrier system
	Do not reuse
	Caution
	Batch code

	Keep dry
	Keep away from sunlight
	Storage temperature limitation
	Do not re-sterilise
	Use-by date
	Do not use if package is damaged and consult the instructions for use
	Catalogue number
	Manufacturer
	Date of manufacture
	Consult the instructions for use
	Carcinogenic, mutagenic, toxic for reproduction or endocrine disruptor
	Medical device
	UDI
	Patient identification
	Name and address of the implanting healthcare institution and provider
	Date of implantation
	Information website for patients

Plugs for baseplate/condyle	100% UHMWPE (ISO 5834-2 and ASTM F 648)
Insert (except UNI full PE), patella	100% UHMWPE (ISO 5834-2 and ASTM F 648)
Full PE insert tibial	99.57 – 99.85 % UHMWPE (ISO 5834-2 and ASTM F 648) + 0.15 – 0.43 % Ti6Al4V ELI (ISO 5832-3 and ASTM F 136)

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