HEXAGON

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FR

NI.D3009A1-2025/01/24

EN

Instructions for use

IT

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ES

Nota de instrucciones

PT

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0. REVISION

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1. DESCRIPTION

> The HEXAGON knee joint prostheses are implantable artificial devices for a knee joint, intended to replace all the articulating surfaces of the damaged/degenerative knee joint (primary total knee replacement) in skeletally mature patients, to reduce pain, and improve quality of life and knee mobility. They shall be implanted by healthcare professional.

> The HEXAGON knee prostheses include femoral, tibial, and patellar components that articulate.

> HEXAGON prostheses are available in four configurations:

- Standard Fixed bearing configuration
- Standard Mobile bearing configuration
- Posterior-Stabilized Fixed bearing configuration
- Posterior-Stabilized Mobile-bearing configuration

Posterior stabilized versions are for use without the posterior cruciate ligament when additional stability is required to prevent anterior subluxation of the femur relative to the tibia in flexion.

Mobile versions allow tibial internal/external rotation during flexion.

These versions are available in cemented or cementless versions.

> The cemented components are intended to be used with bone cement.

> Resurfacing of the patella is not mandatory. Patella is intended to be cemented.

> An optional extension keel can be used to extend the tibial baseplate length.

> The materials used to manufacture these products satisfy current standard requirements (see the end of this leaflet).

2. INDICATIONS

The indications for total knee arthroplasty are as follows:

- Degenerative arthritis, rheumatoid arthritis or

post-traumatic arthritis;

- Post-traumatic loss of the configuration and function of the knee joint;
- Avascular necrosis of the knee;
- Moderate varus, valgus, or flexion deformity in which the ligamentous structures can be returned to adequate function and stability;
- Failure of previous surgical treatment.

3. CONTRAINDICATIONS

Knee replacement is contraindicated in the following cases:

- Any active or suspected latent infection inside or around the knee joint;
- Any neuromuscular or mental disorder which would create an unacceptable risk of instability of the prosthesis, failure of the prosthesis fixation or complications during postoperative care;
- Bone stock compromised due to a disease, an infection or prior implantation which may not provide an adequate support and/or fixation for the prosthesis;
- Skeletal immaturity;
- Severe instability of the knee joint secondary to the lack of integrity and function of the collateral ligaments;
- Hypersensitivities or allergies to the materials used.

4. WARNINGS

> The use of this product is reserved for surgeons with sufficient skills and experience in knee replacement surgery. The product must be used under the conditions specified by the manufacturer for the designated purpose.

> Prior to use, surgeons shall be trained to the assembly and installation techniques and shall have read EUROS documentation. Insertion must be carried out in accordance with the recognised surgical technique. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly.

> Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and the selection and proper positioning of devices, are important factors in the successful

use of this device or system by the surgeon. Furthermore, the selection of a suitable patient together with the latter's cooperation are essential factors for a successful outcome.

Patients who present contra-indications must be turned down.

> The surgeon should take into consideration factors that may affect the performance of the device or system:

- The patient's weight, activity level and occupation ;
- The patient's bone state (osteoporosis, tumors) ;
- The patient's allergies ;
- The patient's diseases and disorders (infectious diseases, mental or neuromuscular disorders) ;
- The patient's addiction (drug, alcoholism or nicotine).

The surgeon must assess the patient's physical and mental state and decide whether the latter is suited to having the device fitted, can stand up to the operation and follow a suitable rehabilitation program.

> The manufacturer is not liable for any complications arising from incorrect diagnosis, incorrect choice of device, the incorrect assembly of components of the device and/or incorrect operating techniques, the limitations of treatment methods or inadequate aseptic conditions.

> Reporting an incident:

Any healthcare professional wishing to make a complaint or express their dissatisfaction in terms of the quality, identification, reliability, safety, efficacy and/or the performance of EUROS device must inform EUROS.

In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email (complaints@euros.fr). The competent authority of the Member State in which the user/or patient is established must be notified in case of a serious incident in relation to the device.

For any additional information involving the devices, please contact EUROS or your distributor.

EUROS SAS

Z.E. Athélia III - 13600 La Ciotat - France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

www.euros.fr

5. PRECAUTIONS

> Caution should be exercised when handling the devices. The devices must not be scratched or damaged in any way whatsoever. Damaged devices and those with damaged packaging must not be used and should be returned to EUROS.

> Devices must be kept in its original sealed packaging in a clean, dry location. Store the sterile products between 10 and 30°C. During transport, temperatures of up to 40°C for short periods can be tolerated. The packaging must not be directly exposed to the sun's rays or extreme temperatures.

> Devices are for single use only and must not be reused. The reuse of a single use device that has come into contact with blood, bone, tissue or other body fluids may cause injury to the patient or user. The risks associated with the reuse of a single-use device may include, but are not limited to, mechanical failure and the transmission of infectious agents.

> Devices must only be used with the instruments specified by EUROS in the reference surgical technique.

> In the event of any disposal due to an error of storage or use, or to excision, devices should be disposed of following the disposal method for hospital waste in accordance with the procedures in force in the establishment.

> Using all means available such as X-rays and templates, the surgeon must ensure that, by guaranteeing the concordance of scales between the X-rays and the templates, the implantation of the device will cause no risks of impingement with other implants already in the patient's body.

Tibial baseplate

> The tibial intramedullary orientation can only be achieved if the tibial curvature allows it.

> The correct selection as well as the correct seating/placement of the prosthetic implant is extremely important for a successful result.

> The posterior slope of the cutting plane must be equal to or higher than 0°. The tibial cutting guide includes a 5° posterior slope.

> A saw blade 1.27 mm thick must be used.

Tibial insert

The tibial insert must correspond to the size of the femoral implant (see the implant matching table).

Femoral condyle

- > Planning is strongly recommended in order to determine the physiological valgus.
- > The implant size must match the size of the anatomical condyle.
- > The HEXAGON ranges are composed of right and left condyles. The sides must be respected.
- > It is important to center the femoral prosthetic implant medio-laterally to avoid any risk of impingement with the lateral ligaments.

Patella

Measure the thickness of the patella to be prosthetized, since regardless of the type of patella (buried or resurfacing), the remaining bone stock after preparation should not be less than 10 mm.

Resurfacing patella

- > Make the bone cut parallel with the patella sagittal axis.
- > The implant must match the diameter of the patella to be prosthetized.

6. EXPECTED PERFORMANCE AND BENEFITS

> Technical performances

- Enabling kinematics and knee joint mobility:
 - Restoring normal physiological ligament tension (lateral and/or posterior cruciate ligament),
 - Being adapted to the different anatomies thanks to a large size range,
 - Restoring the internal rotation of the tibia during movement (extension / flexion),
 - Restoring the dynamic variation of the varus valgus angle during movement (extension / flexion).
- Enabling joint stability:
 - Primary stability: stability of bone anchors,
 - Secondary stability: solidity and durability of cement or coating.
- Resisting to mechanical loads under normal conditions of use:
 - Compression resistance,
 - Integrity of connections between components,
 - Subluxation/dislocation.

- Resisting to wear under normal conditions of use;
- Being adapted to the bony parts ;
- Answering the implantation constraints:
 - Biocompatible,
 - Being easily cleaned and sterilized,
 - Being easily implantable and removable.

> Clinical performances

- Improving knee mobility (depending on the indication) ;
- Reducing pain level during daily activities.

> Clinical benefits

- Decreased pain ;
- Increased quality of life ;
- Easier walking.

7. ADVERSE EFFECTS / RESIDUAL RISKS

The list of potential adverse events associated with any knee arthroplasty includes :

- Loosening or breakage/damage of knee implant components ;
- Reaction of/damage to the tissues, tendons, ligaments or muscles surrounding the prosthesis, blood loss ;
- Dislocation and/or instability of the joint or the prosthesis;
- Incorrect alignment or malposition of implant components (patellar conflict, hypercorrection, misalignment...) ;
- Implant subsidence ;
- Bone fracture or resorption ;
- Nervous disorders ;
- Superficial wound issues, hematoma, swelling, effusion, infection or inflammation ;
- Leg length discrepancy ;
- Low range of motion, stiffness ;
- Pain ;
- Excessive implant wear, bone osteolysis or geode ;
- Heterotopic ossification;
- Cardiovascular disorders, venous thromboembolic disease or pulmonary embolism;
- Patellar problems ;
- Allergy or sensitivity to material ;
- *Intoxication, cellular damage, convulsion ;
- *Cancer, tumor, immunological, and neurological disorders, genetic defects, non-heritable adverse

effects in the progeny and/or an impairment of male or female reproductive functions or capacity.

Treatment of certain undesirable side effects could necessitate additional surgery.

**Note: these residual risks are related to the device biocompatibility. The biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable for the EUROS products.*

8. INFORMATION TO BE PROVIDED TO THE PATIENT

- > The surgeon should inform the patient that the device's safety and durability will be affected by her/his weight, behaviour, namely the physical activities. Before leaving the care facility, the patient must receive instructions and warnings from the surgeon regarding the treatment and therapy to follow, and any restriction on physical and sporting activity or practices. The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, but also of the need for regular medical follow-up. The lifespan of the implant may be shorter than the patient's lifetime or any particular given length of time. The surgeon must inform the patient that the device cannot and does not reproduce the flexibility, strength, reliability or durability of normal healthy bone, and that the device may need to be replaced in the future.
- > Implant card to be provided to the patient: fill the implant card using an indelible ink pen following the mentioned instruction and give it to the patient. Explain him which information are available on the implant card. The surgeon must indicate to the patient that the additional information on the implant is available on the website (accessible with the link shown on the implant card) by means of the product reference, in particular as regards to the precautions for use, the materials or MRI compatibility. The surgeon must inform the patient that they should always have their implant card with them.
- > MR environment: the surgeon must inform the patient to indicate that he is wearing an implant before an MRI examination.

An MRI safety evaluation based on a bibliographical analysis and a demonstration of equivalence has demonstrated the EUROS knee prostheses are MR-conditional. A patient with these devices can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 3.0-Tesla (3.0T) or less ;
- Maximum spatial gradient field of 750-Gauss/cm ;
- Whole-body-averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning for patient in a 3T MR system.

MR imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the EUROS knee device.

- > The surgeon must inform the patient to report any incident in relation to the devices to the surgeon.
- > CMR substances: the surgeon must inform the patient that the HEXAGON condyles and baseplates contain more than 0.1% (weight/weight) of cobalt identified as a carcinogenic and toxic for reproduction substance according to the regulation 1272/2008, 1907/2006 and 2020/2096. However, the cobalt chromium alloy has long clinical history and remains the best alternative in knee prosthesis. Our biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable.
- > The surgeon should inform the patient about the residual risks available in section 7 of this instruction for use.

9. CHOICE OF DEVICE / DEVICE COMPATIBILITY

- > The surgeon is responsible for selecting the devices according to the patient's morphology and the nature of the condition.
This choice may require preoperative planning. The instrumentation provided with the implant can also guide this choice thanks to a range of trial implants corresponding to the range of implants.
- > Information on device size is provided on the device packaging and on the device itself. The surgical team must check this information before the surgery.
- > The use of a cemented or cementless implant depends on the surgeon's choice.

List of cemented/cementless implants of products:

	Cemented	Cementless
Condyle	✓	✓
Baseplate	✓	✓
Patella	✓	
Tibial extension keel		✓

> In the case of devices designed to be implanted with cement or bone substitute, no incompatibility of use has been detected. The selection of the cement is the discretion of the surgeon.

Regarding devices designed to be implanted without cement, the manufacturer cannot, under any circumstances, be held responsible for implanting the device with cement.

COMBINATION OF DEVICES

> EUROS devices must be combined only with EUROS devices. All combination with devices of other manufacturer is forbidden.

> Compatibility tables representing the authorised combinations between the components for correct assembly:

HEXAGON FIXE

Fixed baseplate	Fixed insert	Femoral component sizes									
		0	1	2	3	4	5	6	7	8	9
0	0	✓	✓	✓							
1	1/2	✓	✓	✓	✓	✓					
2											
3	3/4		✓	✓	✓	✓	✓	✓			
4											
5	5/6				✓	✓	✓	✓	✓	✓	
6											
7	7/8						✓	✓	✓	✓	✓
8											

HEXAGON mobile

Femoral component sizes	Mobile insert	Mobile baseplate sizes									
		0	1	2	3	4	5	6	7	8	
0	0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
1	1	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
2	2		✓	✓	✓	✓	✓	✓	✓	✓	✓
3	3			✓	✓	✓	✓	✓	✓	✓	✓
4	4				✓	✓	✓	✓	✓	✓	✓
5	5					✓	✓	✓	✓	✓	✓
6	6						✓	✓	✓	✓	✓
7	7							✓	✓	✓	✓
8	8								✓	✓	✓
9	9									✓	✓

10. STERILITY

> EUROS devices are supplied sterile (gamma irradiation) at a sterility assurance level (SAL) of 10^{-6} .

> Inspect every package before use. Do not use the component if the integrity of the double envelope ensuring the maintenance of sterility has been compromised (loss of vacuum, damaged envelopes, damaged seal area), or if the expiry date has elapsed.

> Once the package has been opened, the device must never be re-sterilised. Any re-sterilisation may have adverse effects on the device materials and cause premature defects.

> For instruments used to insert implants, see the relevant leaflet for information on sterilisation.

11. ADDITIONAL INFORMATION

HEXAGON basic UDI-DI : 37004372G00496

12. MEANING OF THE SYMBOLS

List of symbols applicable for the devices covered by this instruction for use:

	Sterilized using irradiation Double sterile barrier system
	Do not reuse
	Caution
	Batch code
	Keep dry
	Keep away from sunlight
	Storage temperature limitation
	Do not re-sterilise
	Use-by date
	Do not use if package is damaged and consult the instructions for use

	Catalogue number
	Manufacturer
	Date of manufacture
	Consult the instructions for use
	Carcinogenic, mutagenic, toxic for reproduction or endocrine disruptor
	Medical device
	UDI
	Patient identification
	Name and address of the implanting healthcare institution and provider
	Date of implantation
	Information website for patients

13. MEANING OF THE ABBREVIATIONS

Qty: quantity - fr: French - en: English - it: Italian - es: Spanish

CoCr: Cobalt chromium - Ti: titanium - UHMWPE : ultra high molecular weight polyethylene

14. MATERIAL

Cementless condyle	96.65% – 96.74% CoCr (ISO 5832-4 and ASTM F75) + 3.26% – 3.35% Ti coating (ISO 13179-1 and ASTM F1580)
Cementless baseplate	93.21% – 95.98% CoCr (ISO 5832-4 and ASTM F75) + 0.09% – 0.1% Ti coating (ISO 13179- and ASTM F1580) + 3.91%– 6.70% Ti6Al4V ELI (ISO 5832-3 and ASTM F136)
Cemented condyle	100% CoCr (ISO 5832-4 and ASTM F75)
Cemented baseplate	90.64% – 94.97% CoCr (ISO 5832-4 and ASTM F75) + 5.03% – 9.36% Ti6Al4V ELI (ISO 5832-3 and ASTM F136)
Insert, patella	100% UHMWPE (ISO 5834-2 and ASTM F648)
Keel	100% Ti6Al4V ELI (ISO 5832-3 and ASTM F136)

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