SCULTRA II

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FR

NI.D4002A1-2024/07/11

EN

Instructions for use

IT

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ES

Nota de instrucciones

DE

Gebrauchsanweinsug

NL

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PT

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0. REVISION

Version & date	Modifications	URL internet
A1 - 11/07/2024	Creation,	-

1. DESCRIPTION

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1. DESCRIPTION

- > The SCULTRA II shoulder prostheses are intended to reconstitute the glenohumeral joint in skeletally mature patients with humeral fracture, severe shoulder arthropathy or a previously failed shoulder joint replacement.
- > The range consists of the following implants :
 - SCULTRA II STANDARD (Anatomical shoulder prosthesis) composed of a humeral stem, offset ring, angulation cones, humeral head and glenoid;
 - SCULTRA II REVERSE (Reverse shoulder prosthesis) including the REVISION variant. SCULTRA II REVERSE is composed of a stem, insert, glenoid baseplate, glenoid sphere. A raising block may also be used.
 - SCULTRA II REVERSE REVISION is composed of a stem, hooded humeral head and spacer with locking screws;
 - Cancellous Screw Ø 4.5 mm, Cancellous Screw Ø 6.5 mm.
- > The SCULTRA II STANDARD variant shoulder prosthesis can be used in total or hemi shoulder arthroplasty.
- > The SCULTRA II REVERSE variant shoulder prosthesis is used in total reverse shoulder arthroplasty combined with a grossly rotator cuff deficient shoulder joint.
- > The SCULTRA II REVERSE REVISION variant is intended to be used in hemi cuff-tear shoulder arthroplasty.
Hemi shoulder joint prostheses are only intended to be used in cases where the natural glenoid of the patient is in good condition.
- > The materials used to manufacture these products satisfy current standard requirements (see the end of this leaflet).

2. INDICATIONS

SCULTRA II may be indicated to relieve severe pain or significant disability caused by:

- Degenerative arthritis (i.e. osteoarthritis), rheumatoid arthritis, or traumatic arthritis of the glenohumeral joint,
- Traumatic disease of the glenohumeral joint,
- Non-union humeral head fracture,
- Avascular necrosis of the humeral head,
- Fracture-dislocations of the proximal humerus where the articular surface is severely comminuted, separated from its blood supply or where the surgeon's experience indicates that alternative methods of treatment are unsatisfactory,
- Other difficult clinical problems where shoulder arthrodesis or resection arthroplasty are not acceptable,
- Revision surgery where other treatments or devices have failed.

The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the subject device.

SCULTRA II REVERSE

SCULTRA II REVERSE is indicated for patients with a grossly rotator cuff deficient shoulder joint with severe arthropathy or a previously failed joint replacement with a grossly rotator cuff deficient shoulder joint.

3. CONTRAINDICATIONS

Shoulder joint replacement is contraindicated in the following cases:

- Absence or severe impairment of deltoid contraction,
- Recent or active local or systemic infection, sepsis and osteomyelitis,
- Conditions that would eliminate or tend to eliminate adequate implant and/or cement support such as insufficient quantity or quality of bone (particularly glenoid) resulting from conditions such as osteoporosis, trauma, or neoplasia,
- Neuropathic joint,
- Metabolic disorders which may impair bone formation and/or wound healing,
- Allergy to one of the chemical compounds,
- Signs of local inflammation,
- Absence of covering tissues,
- Pathological obesity,

- Mental pathologies,
- Implant overall dimensions incompatible with the size of the receiving site,
- Allergy or intolerance to titanium and to its alloy,
- The age and the physical tests of the patient are not compatible with the surgery.

SCULTRA II STANDARD

SCULTRA II STANDARD is contraindicated in the cases of irreparable rotator cuff tear, and cases where stabilizing musculature cannot provide stability.

4. WARNINGS / CAUTIONS

- > The use of this product is reserved for surgeons with sufficient skills and experience in shoulder replacement surgery. The product must be used under the conditions specified by the manufacturer for the designated purpose.
- > Prior to use, surgeons shall be trained to the assembly and installation techniques and shall have read EUROS documentation. Insertion must be carried out in accordance with the recognised surgical technique. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly.
- > Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and the selection and proper positioning of devices, are important factors in the successful use of this device or system by the surgeon. Furthermore, the selection of a suitable patient together with the latter's cooperation are essential factors for a successful outcome.
Patients who present contra-indications must be turned down.
- > The surgeon should take into consideration factors that may affect the performance of the device or system:
 - The patient's weight, activity level and occupation
 - The patient's bone state (osteoporosis, tumors)
 - The patient's allergies,
 - The patient's diseases and disorders (infectious diseases, mental or neuromuscular disorders)
 - The patient's addiction (drug, alcoholism or nicotine)

The surgeon must assess the patient's physical and mental state and decide whether the latter is suited to having the device fitted, can stand up to the

operation and follow a suitable rehabilitation program.

- > The manufacturer is not liable for any complications arising from incorrect diagnosis, incorrect choice of device, the incorrect assembly of components of the device and/or incorrect operating techniques, the limitations of treatment methods or inadequate aseptic conditions.

- > Reporting an incident:

Any healthcare professional wishing to make a complaint or express their dissatisfaction in terms of the quality, identification, reliability, safety, efficacy and/or the performance of EUROS device must inform EUROS.

In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email (complaints@euros.fr). The competent authority of the Member State in which the user/or patient is established must be notified in case of a serious incident in relation to the device.

For any additional information involving the devices, please contact EUROS or your distributor.

EUROS SAS

Z.E. Athélia III - 13600 La Ciotat - France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

www.euros.fr

5. PRECAUTIONS

- > Caution should be exercised when handling the devices. The devices must not be scratched or damaged in any way whatsoever. Damaged devices and those with damaged packaging must not be used and should be returned to EUROS.
- > Devices must be kept in its original sealed packaging in a clean, dry location. Store the sterile products between 10 and 30°C. During transport, temperatures of up to 40°C for short periods can be tolerated. The packaging must not be directly exposed to the sun's rays or extreme temperatures.
- > Devices are for single use only and must not be reused. The reuse of a single use device that has come into contact with blood, bone, tissue or other body fluids may cause injury to the patient or user. The risks associated with the reuse of a single-use device may include, but are not limited to,

mechanical failure and the transmission of infectious agents.

- > Devices must only be used with the instruments specified by EUROS in the reference surgical technique.
- > In the event of any disposal due to an error of storage or use, or to excision, devices should be disposed of following the disposal method for hospital waste in accordance with the procedures in force in the establishment.
- > Using all means available such as X-rays and templates, the surgeon must ensure that, by guaranteeing the concordance of scales between the X-rays and the templates, the implantation of the device will cause no risks of impingement with other implants already in the patient's body.

SCULTRA II STANDARD

- > For most cases, glenoids of size Small should be used with all the humeral heads. Medium and large glenoids should be reserved to extreme cases.
- > A radio-opaque marker has been inserted into the glenoid implant to enable radiographic verification of the correct positioning.
- > Prior to implanting the humeral head, it is important to set the ring to the selected center off-set degree (See SCULTRA II STANDARD Surgical Technique).

SCULTRA II REVERSE

- > It is imperative to respect the correct positioning of the base so as to enable the glenoid sphere to be center off-set downwards. A locating pin has been provided for this purpose on the upper pole of the base (See SCULTRA II REVERSE Surgical Technique).

6. EXPECTED PERFORMANCE AND BENEFITS

- > Technical performances
 - Enabling kinematics and shoulder joint mobility:
 - Restoring normal physiological joint tension;
 - Being adapted to the different anatomies thanks to a large size range;
 - Restoring the range of motion;
 - Restoring the appropriate retroversion and neck angulation;
 - Maintaining the center of rotation;
 - Enabling joint stability:

- Primary stability: stability of bone anchors;
- Secondary stability : Solidity and durability of cement or coating;
- Resisting to mechanical loads under normal conditions of use:
 - Locking mechanism resistance;
 - Integrity of connections between components;
 - Subluxation/dislocation ;
- Resisting to wear under normal conditions of use;
- Answering the implantation constraints:
 - Biocompatible;
 - Being easily cleaned and sterilized;
 - Being easily implantable and removable.
- > Clinical performances
 - Restoration of shoulder articulation ;
 - Restoration of shoulder range of motion.
- > Clinical benefits
 - Improvement of patient functional outcome;
 - Pain relief ;
 - Safe and reliable fixation of the prosthesis.

7. ADVERSE EFFECTS / RESIDUAL RISKS

The list of potential adverse events includes :

- Excessive wear of the implant components secondary to impingement of components or damage of articular surfaces;
- Instability;
- Fracture, migration, loosening, dissociation, impingement, subluxation, or dislocation of the prosthesis or any of its components;
- Possible detachment of the coating(s), potentially leading to increased debris particles;
- Bone fracture;
- Scapular notching only for reverse shoulder implants;
- Glenoid erosion only for hemi arthroplasty;
- Osteolysis;
- Soft or surrounding tissues injuries;
- Metal sensitivity reactions or other allergic/histological reactions to implant materials; inflammation, formation of erythema or oedema or implant rejection;
- Superficial or deep infection, septic loosening;

- Delayed wound healing;
- Vascular damage resulting in blood loss and/or hematoma;
- Neurologic injury or neuropathy resulting in transient or permanent dysesthesia, weakness, pain, and/or numbness;
- Muscle or soft tissues injury, secondary rotatory cuff tear;
- Cardiovascular disorders including venous thrombosis, pulmonary embolism, myocardial infarction, cerebrovascular accident or neurovascular complications;
- Decreased range of motion;
- Stiff or frozen shoulder;
- *Intoxication, cellular damage, convulsion
- *Cancer, tumor, immunological, and neurological disorders, genetic defects, non-heritable adverse effects in the progeny and/or an impairment of male or female reproductive functions or capacity

Treatment of certain undesirable side effects could necessitate additional surgery.

**Note: these residual risks are related to the device biocompatibility. The biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable for the EUROS products.*

8. INFORMATION TO BE PROVIDED TO THE PATIENT

> The surgeon should inform the patient that the device's safety and durability will be affected by her/his weight, behaviour, namely the physical activities. Before leaving the care facility, the patient must receive instructions and warnings from the surgeon regarding the treatment and therapy to follow, and any restriction on physical and sporting activity or practices. The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, but also of the need for regular medical follow-up. The lifespan of the implant may be shorter than the patient's lifetime or any particular given length of time. The surgeon must inform the patient that the device cannot and does not reproduce the flexibility, strength, reliability or durability of normal healthy bone, and that the device may need to be replaced in the future.

> Implant card to be provide to the patient: fill the implant card using an indelible ink pen following the mentioned instruction and give it to the patient. Explain him which information are available on the implant card. The surgeon must indicate to the patient that the additional information on the implant is available on the website (accessible with the link shown on the implant card) by means of the product reference, in particular as regards to the precautions for use, the materials or MRI compatibility. The surgeon must inform the patient that they should always have their implant card with them.

> MR environment: the surgeon must inform the patient to indicate that he is wearing an implant before an MRI examination.

An MRI safety evaluation based on a bibliographical analysis and a demonstration of equivalence has demonstrated the EUROS shoulder prostheses are MR-conditional. A patient with these devices can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0-Tesla
- Maximum spatial gradient field of 720-Gauss/cm
- Whole-body-averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning per pulse sequence.

MR imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the EUROS shoulder prosthesis.

> The surgeon must inform the patient to report any incident in relation to the devices to the surgeon.
> The surgeon should inform the patient about the residual risks available in section 7 of this instruction for use.

9. CHOICE OF DEVICE / DEVICE COMPATIBILITY

> The surgeon is responsible for selecting the devices according to the patient's morphology and the nature of the condition.

This choice may require preoperative planning. The instrumentation provided with the implant can also guide this choice thanks to a range of trial implants corresponding to the range of implants.

- > Information on device size is provided on the device packaging and on the device itself. The surgical team must check this information before the surgery.
- > The use of a cemented or cementless implant depends on the surgeon's choice.
- > List of cemented/cementless implants of products:

	Cemented	Cementless
SCULTRA II - STANDARD Humeral stems	✓	✓
SCULTRA II - STANDARD Glenoid	✓	
SCULTRA II - REVERSE Humeral stems	✓	✓
SCULTRA II - REVERSE Glenoid baseplate		✓

- > In the case of devices designed to be implanted with cement or bone substitute, no incompatibility of use has been detected. The selection of the cement is the discretion of the surgeon.
Regarding devices designed to be implanted without cement, the manufacturer cannot, under any circumstances, be held responsible for implanting the device with cement.

COMBINATION OF DEVICES

- > EUROS devices must be combined only with EUROS devices. All combination with devices of other manufacturer is forbidden.
- > All components SCULTRA II STANDARD (of any size) are compatible with one another.
- > All components SCULTRA II REVERSE (of any size) are compatible with one another.
- > All components SCULTRA II REVISION (of any size) are compatible with one another.

10. STERILITY

- > EUROS devices are supplied sterile (gamma irradiation) at a sterility assurance level (SAL) of 10^{-6} .
- > Inspect every package before use. Do not use the component if the integrity of the double envelope ensuring the maintenance of sterility has been compromised (loss of vacuum, damaged envelopes, damaged seal area), or if the expiry date has elapsed.
- > Once the package has been opened, the device must never be re-sterilised. Any re-sterilisation may have

adverse effects on the device materials and cause premature defects.

- > For instruments used to insert implants, see the relevant leaflet for information on sterilisation.

11. ADDITIONAL INFORMATION

The Summary of Safety and Clinical Performance is available upon request (e-mail: qualite@euros.fr ; Tel : +33 4 42 71 42 71).

SCULTRA II basic UDI-DI : 37004372E0018J

12. MEANING OF THE SYMBOLS

List of symbols applicable for the implants covered by this instruction for use

	Sterilized using irradiation Double sterile barrier system
	Do not reuse
	Caution
	Batch code
	Keep dry
	Keep away from sunlight
	Storage temperature limitation
	Do not re-sterilise
	Use-by date
	Do not use if package is damaged and consult the instructions for use
	Catalogue number
	Manufacturer
	Date of manufacture

	Consult the instructions for use
	Medical device
	UDI
	Patient identification
	Name and address of the implanting healthcare institution and provider
	Date of implantation
	Information website for patients

13. MEANING OF THE ABBREVIATIONS

Qty: quantity - fr: French - en: English - it: Italian - es: Spanish - de: German - nl: Dutch

HAp: Hydroxyapatite powder - Ti6Al4V ELI : aluminum-vanadium alloy titanium - UHMWPE: ultra high molecular weight polyethylene - X2CrNiMo 18-15-3 : Stainless steel 316L

14. MATERIAL

SCULTRA II STANDARD

SCULTRA II STD - Humeral stem - cemented	100% Ti6Al4V ELI (ISO 5832-3 and ASTM F620)
SCULTRA II STD - Humeral stem - cementless	97.58 - 98.92% Ti6Al4V ELI (ISO 5832-3 and ASTM F620) + 1.08 - 2.42% HAp (ISO 13779-6)
SCULTRA II STD - Cone	100% Ti6Al4V ELI (ISO 5832-3 and ASTM F136)
SCULTRA II STD - Offset ring	100% Ti6Al4V ELI (ISO 5832-3 and ASTM F136)
SCULTRA II STD - Humeral head	100% X2CrNiMo 18-15-3 (ISO 5832-1 and ASTM F138)
SCULTRA II STD - Glenoid cemented and radiopaque marker	99.78% UHMWPE (ISO 5834-2 and ASTM F648) + 0.22% X2CrNiMo 18-15-3 (ISO 5832-1 and ASTM F138)

SCULTRA II REVERSE

SCULTRA II Reverse - Humeral stem - cemented	100% Ti6Al4V ELI (ISO 5832-3 and ASTM F620)
--	---

SCULTRA II Reverse - Humeral stem - cementless	96.47 - 98.45% Ti6Al4V ELI (ISO 5832-3 and ASTM F620) + 1.55 - 3.53% HAp (ISO 13779-6)
SCULTRA II Reverse - Spacer with locking screw	100% Ti6Al4V ELI (ISO 5832-3 and ASTM F136)
SCULTRA II Reverse - Raising block with locking screw	100% Ti6Al4V ELI (ISO 5832-3 and ASTM F136)
SCULTRA II Reverse - Glenoid Sphere with Locking Screw	93.61% X2CrNiMo 18-15-3 (ISO 5832-1 and ASTM F138) + 6.39% Ti6Al4V ELI (ISO 5832-3 and ASTM F136)
SCULTRA II Reverse - Glenoid baseplate cementless	95.36% Ti6Al4V ELI (ISO 5832-3 and ASTM F136) + 4.64% HAp (ISO 13779-6)
SCULTRA II Reverse - Humeral insert	100% UHMWPE (ISO 5834-2 and ASTM F648)
SCULTRA II Reverse - Revision Hooded humeral head	100% X2CrNiMo 18-15-3 (ISO 5832-1 and ASTM F138)
SCULTRA II Reverse - Revision glenoid baseplate	97.97% X2CrNiMo 18-15-3 (ISO 5832-1 and ASTM F138) + 2.03% Ti6Al4V ELI (ISO 5832-3 and ASTM F136)

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