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**E.SPINE® - TORQUE WRENCH - B2524052**

FOR MULTIPLE USE - NOT STERILE


INSTRUCTIONS FOR USE*FOR SURGEONS AND MEDICAL STAFF***1 / PURPOSE**

The E.SPINE® torque wrench B2524052 is an instrument to tighten the locking screws on the E.SPINE® posterior stabilization system to 10Nm ($\pm 10\%$).

2 / DESCRIPTION

The E.SPINE® torque wrench B2524052 is comprised of the following parts:

- An ergonomic handle for a firm grip and to ensure the limit torque is reached
- A system that releases at a torque limit of 10Nm
- A ¼" female connector for assembly onto a ¼" male

The metal parts are made from stainless steel. This instrument effectively locks the E.SPINE® posterior stabilization system containing metal rods.

3 / WARNINGS / CAUTIONS

It is strongly recommended that the device be implanted only by surgeons familiar with the general problems of spinal surgery and who are proficient in the surgical techniques specific to the product. Lack of experience or expertise with this instrument can cause complications.

It is the surgeon's responsibility to familiarise themselves with the assembly and installation techniques prior to use, by obtaining the documentation available from EUROS representatives and distributors. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly

The manufacturer is not liable for any complications due to improper use of the instrument, incorrect assembly and/or incorrect operating techniques or inadequate sterilization.

4 / PRECAUTIONS FOR USE

Handle the instrument with care. This instrument should only be used for locking appropriate implants. Any other use may damage the instrument.

Check that all parts are present and in good condition before use.

Any impact to the instrument may cause a malfunction (fall, use of a hammer or impaction instrument on the handle, etc.). Handles that have received an impact must be returned to EUROS.


Handles damaged in any other way must not be used and must be returned to EUROS.

This instrument should only be used with the instruments and implants specified by EUROS for the reference surgical technique and in the «Implant/instrument compatibility» section of these instructions.

5 / IMPLANT/INSTRUMENT COMPATIBILITY

The E.SPINE® torque wrench B2524052 must be used with:

- E.SPINE® locking screwdriver B2522031
- E.SPINE® counter torque wrench B2522051

The torque wrench used with the locking screwdriver and counter torque wrench must be used only for tightening the following locking screws on the E.SPINE® posterior stabilization system:

E.SPINE®
Hexagonal locking screw - 4.5 mm
Reference A2511001

The E.SPINE® torque wrench B2524052 can only be used with compatible instruments recommended by the manufacturer and must not be used with other instruments with a ¼" connector.

6 / CLEANING, DECONTAMINATION AND STERILIZATION

All E.SPINE® instruments are delivered non-sterile. It is the responsibility of the healthcare facility to clean, decontaminate and sterilize the instrument prior to use, in accordance with approved methods.

If the device is delivered in packaging, the plastic envelope must be removed before any procedure.

Decontamination and cleaning must be carried out before sterilization. Decontamination reduces the number of microorganisms and facilitates subsequent cleaning. After cleaning and decontamination, the instrument must be rinsed carefully with water.

In order to maintain optimal functionality and safety of the device, the following must be avoided: use of metal brushes, scouring pads or objects that could impair the device, use of chlorine-based chemicals, acids or organic or ammonia-based solvents. Sodium hypochlorite should not be used as it corrodes stainless steel.

Washer/disinfecter cycle:

It is recommended that a cycle be run that is specifically for the inactivation of prions.

| | |
|-----------------------|---|
| 1 st phase | Pre-wash: softened cold water - 2min |
| 2 nd phase | Washing – prion inactivation: softened hot water, washing with a detergent, t°55-65°C - 10min |
| 3 rd phase | Neutralisation: softened hot water, neutralising agent - 3min |
| 4 th phase | First rinse: softened water - 3min |
| 5 th phase | Second rinse: RO water - 2min drying activator, Thermal disinfection: t°93°C - 3min |
| 6 th phase | Drying |

In the case of patients with suspect or diagnosed TSE, the cleaning and decontamination procedure must be carried out in accordance with the instructions of Directive DGS/RI3/2011/449.

Sterilization cycle for medical devices:

Local or national specifications must be followed where steam sterilization requirements are stricter or more conservative than those listed in the table.

Always sterilize surgical instruments using a container suitable for steam sterilization and follow the instructions below:

| Method | Type of Cycle | Pressure | Temperature | Exposure time | Minimum drying time ² |
|--------------------|---------------|----------|-------------|---------------|----------------------------------|
| Steam ¹ | Pre-vacuum | 3 bars | 134°C | 18 min | 30 min |

¹ Steam sterilization parameters recommended in accordance with Directive DGS/RI3/2011/449.

² Minimum drying times have been validated using sterilizers equipped for vacuum drying. Drying cycles using ambient atmospheric pressure may require longer drying times. Please refer to the manufacturer's recommendations.

NOTE :

The sterilizer manufacturer's operating instructions and instructions for load configuration must be followed absolutely.

7 / CARE AND MAINTENANCE OF THE INSTRUMENT

Follow the maintenance instructions below for stainless steel surgical instruments:

- Do not allow blood or tissue to dry on the instrument
 - Rinse the instrument immediately after use before decontamination
 - If possible, separate instruments made from different metals
- Check that each instrument is clean and functioning properly before sterilizing
Check the release threshold yearly.

The instrument's lifespan is estimated at 6,000 uses, i.e. 3 years' intensive use of the handle. After this, the instrument is no longer guaranteed to function properly and it must not be used.

8 / ADDITIONAL INFORMATION

Any healthcare professional who has a complaint or is dissatisfied with the quality, identification, reliability, safety, efficiency and/or performance of the E.SPINE® torque wrench B2524052 must notify EUROS.

In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email (qualite@euros.fr).

For any further information about the E.SPINE® device, please contact EUROS or your distributor.

EUROS SAS

Z.E. Athélia III - 13600 La Ciotat - France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

www.euros-orthopaedics.com

9 / HOW TO READ THE SYMBOLS

Warning: check the instruction leaflet and the operating technique

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