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(En) — NLD1009-1010-1021B - 2019/11/26

INSTRUCTIONS FOR USE FOR SURGEONS AND MEDICAL STAFF

E.2S IMPLANTS

1. PURPOSE

E.2S cervical plates are implants designed to be placed using an anterior approach in order to join from 2 to 6 cervical vertebrae together.

2. DESCRIPTION

> The implant is held in place by means of fixing screws, with two being placed in each vertebral body.

The cervical plates are composed of two separate elements:

- The cervical plate itself,
- Corpeal fixing screws.

> The E.2S cervical plates come in different lengths, with a varying number of holes depending on the number of vertebrae to be instrumented. The screws designed for use with the E.2S plates are available in different diameters and lengths.

> The E.2S cervical plates and their associated screws are made of titanium alloy (TA6V4 ELI) in accordance with ISO 5832-3 and ASTM F 136.

3. INDICATIONS

The list (not exhaustive) of indications includes:

- Trauma
- Degenerative disorders of the spine, including disc diseases, spondylolisthesis, hernia, stenosis etc.
- Tumours
- Failure of a previous operation (pseudarthrosis)

These indications must be considered during discussion with the patient. This list is not exhaustive and only the surgeon can make a decision.

4. CONTRAINDICATIONS

The list (not exhaustive) of contraindications includes :

- Signs of local inflammation or infection
- Severe osteoporosis or severe bone resorption
- Fever or leucocytosis
- Pathological obesity
- Pregnancy
- Patients whose physical, mental or medical condition would prevent a positive surgical outcome
- Suspected or proven allergy or intolerance to the implanted materials
- Use in combination with other metal systems
- Patients having inadequate tissue coverage
- Patients who do not wish to restrict their activities or follow medical advice

These contraindications must be considered during discussion with the patient. This list is not exhaustive and only the surgeon can make a decision.

5. WARNINGS/ CAUTIONS

> It is strongly recommended that the device be implanted only by surgeons familiar with the general problems of spinal surgery and who are proficient in the surgical techniques specific to the product. Lack of experience or expertise with these implants can cause complications.

> It is the surgeon's responsibility to familiarise themselves with the assembly and installation techniques prior to use, by obtaining the documentation available from EUROS representatives and distributors. Implanting must be carried out according to the instructions of the recommended surgical procedure. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly.

> Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and the selection and proper positioning of implants, are important factors in the successful use of this implant or system by the surgeon. Furthermore, appropriate patient selection and patient cooperation will significantly affect results. Patients subject to contraindications should be avoided.

> The physician should take into consideration factors that may affect the performance of the implant or system:

- Obesity
- Malnutrition
- Osteoporosis or osteopenia
- Unsuitable professional, sports-related or [other] physical activity
- Bone tumour at implant anchoring sites
- Addiction, chemical dependency, alcoholism or smoking
- Infection, infectious disease
- Signs of local inflammation
- Trauma
- Pregnancy
- Mental illness
- History of falls
- (Allergic) reaction to foreign bodies
- Extremely distorted anatomy due to birth defects

- Fever or leukocytosis
- Insufficient or inadequate bone or muscle quality
- Nerve palsy

> Patients should be informed and warned of the consequences arising from this.

> The surgeon must inform the patient who is to receive the implanted device that the safety and durability of the implant depends on their weight, their behaviour and particularly their physical activity. Before leaving the care facility, the patient must receive instructions and warnings from the surgeon with regard to the treatment and therapy to follow, and to any restriction on their physical and sporting activity or practices. The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, but also of the need for regular medical follow-up monitoring. The lifespan of the implant may be shorter than the patient's lifetime or any particular given length of time. The surgeon must inform the patient that the device cannot and does not reproduce the flexibility, strength, reliability or durability of normal healthy bone, and that the device may need to be replaced in the future.

> The manufacturer is not liable for any complications arising from incorrect diagnosis, incorrect choice of implant, the incorrect assembly of components of the implant and/or incorrect operating techniques, the limitations of treatment methods or inadequate asepis.

NOTE TO PHYSICIAN:

Although it is the doctor who is in possession of the [relevant] knowledge and who is the intermediary between the manufacturer and the patient, the important medical information in this leaflet must be conveyed to the patient.

6. PRECAUTIONS

> Use caution when handling the implants. The implants must not be scratched or damaged in any other way. Damaged implants and those whose packaging is damaged must not be used and should be returned to EUROS.

> Implants should be kept in their original sealed packaging in a clean, dry place. This packaging must not be exposed to the direct action of sunlight, ionizing radiation, extreme temperatures or particle contamination.

> This device is for single use only and should not be reused. The reuse of a single-use device that has come into contact with blood, bone, tissue or other body fluids may cause injury to the patient or user. The risks associated with the reuse of a single-use device may include, but are not limited to, mechanical failure and the transmission of infectious agents.

> Implants must only be used with the instruments specified by EUROS in the reference surgical technique.

> In the event of any disposal due to an error of storage or use, or to excision, implants should be disposed of following the disposal method for hospital waste in accordance with the procedures in force in the establishment.

7. ADVERSE EFFECTS

The list (not exhaustive) of adverse effects includes:

- Migration of the implant
- Breakage of the implant
- Patient's reaction to a foreign body
- Instability and degeneration of adjacent segments
- Pseudarthrosis
- Vertebral fracture
- Superficial and/or deep infection
- Injury to surrounding soft tissues (vascular, nervous, visceral, muscular)
- Dural leakage
- Neurological injury
- Paralysis
- Surgical complications (thrombosis, pulmonary embolism, myocardial infarction, bursitis, seroma, haematoma, difficulty healing)
- Pain
- Death

Some adverse effects may require additional surgery.

8. CHOICE OF IMPLANT / IMPLANT COMPATIBILITY

> The surgeon is responsible for the choice of implants, depending on the patient and the size of their vertebrae. This choice may require preoperative planning. The instrumentation provided with the implant can also guide this choice thanks to a range of trial implants corresponding to the range of implants.

> Information about the size and shape of the implant is shown on the wrapping of the device. The surgical team must check this information before implantation.

9. STERILITY

> The cervical plates are supplied sterile (gamma irradation) to a sterility assurance level (SAL) of 10⁻⁶.

Inspect every package before use.

Do not use the component if the integrity of the double envelope ensuring the maintenance of sterility has been compromised (loss of vacuum, damaged envelopes, damaged seal area), or if has passed its expiry date.

> Once the package has been opened, the implant must never be re-sterilized. Any re-sterilization may have undesirable effects on the materials of the implant and cause premature failure.

> For the instrumentation associated with the installation of the cervicale plates, refer to the corresponding leaflet concerning sterilization information.

10. MRI INFORMATION

The safety and compatibility of the E.2S devices have not been evaluated in a magnetic resonance (MR) environment. The heating, migration or image artefacts for the E.2S device have not been evaluated in a magnetic resonance (MR) environment.

MRI examination of a patient fitted with this device can cause injury to the patient.

11. ADDITIONAL INFORMATION

Any health care professional who has a complaint or cause for dissatisfaction with regard to the quality, identification, reliability, safety, efficiency and/or performance of E.2S devices cages must inform EUROS.

In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email (qualite@euors.fr).

For any further information about the E.2S devices, please contact EUROS or your distributor.

EUROS SAS

Z.E. Athélla III – 13600 La Ciotat – France

Tel : +33 4 42 71 42 71 | Fax : +33 4 42 71 42 80

www.euors-orthopaedics.com

12. HOW TO READ THE SYMBOLS

	<i>Sterilized by irradiation</i>		<i>Use-by date</i>
	<i>Do not reuse</i>		<i>Consult the precautions for use</i>
	<i>Do not re-sterilise</i>		<i>Batch code</i>
	<i>Warning: check the instruction leaflet and the operating technique</i>		<i>Catalogue reference</i>
	<i>Do not use if packaging is damaged.</i>		<i>Manufacturer</i>
	<i>Keep dry</i>		<i>Date of manufacture</i>
	<i>Keep away from sunlight</i>		

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