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## CROWN VOÛTE CRÂNIENNE

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#### DESTINÉE AUX CHIRURGIENS

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### 1. DESCRIPTION

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## CROWN CRANIAL VAULT

### STERILE CUSTOM-MADE IMPLANT

#### INTENDED FOR SURGEONS AND MEDICAL PERSONNEL

#### 1. DESCRIPTION

> The Crown cranial vault is a custom-made device designed and manufactured from the patient's medical anatomical data for which it is intended. Under no circumstances may it be used for a recipient other than the patient for whom it was designed.

> The Crown cranial vault has a constant thickness. It is bevelled on its contour to allow impingement-free positioning with the bone. A hole mesh is present on the entire surface to allow the surgeon to suspend the dura mater and reduce the vacuum spaces between the implant and the corpus callosum.

> Le The Crown device is made of 100% implantable polyetheretherketone (PEEK OPTIMA® LT1) according to ASTM F 2026.

#### 2. INTENDED USE

> The Crown device is intended to reconstruct the missing part of the cranial vault.

#### 3. TARGET PATIENT GROUP

> The Crown cranial vault can be used for adult, adolescent or paediatric patients (with a fully-developed cranial cavity).

#### 4. USERS

> Use of this custom device is limited to surgeons with sufficient skills and experience in the cranioplasty field. The product must be used under the conditions set out by the manufacturer and for the purposes intended by the manufacturer.

#### 5. INDICATIONS

> The Crown cranial vault is indicated for the reconstruction of cranial bone defects in the following cases:

- Trauma
- Tumour resection
- Decompressive craniectomy
- Revisions

These indications must be considered during the examination of the patient. Only the surgeon can make a decision.

#### 6. EXPECTED PERFORMANCE AND BENEFITS

##### > TECHNICAL PERFORMANCE

- Allows reconstruction of a cranial vault bone defect
  - Fills inside of the bone defect
  - Rebuilds the aesthetics of the outer table
  - Ensures the stability of the implant on the vault

##### > CLINICAL PERFORMANCE

- Completes haemostasis and prevents extradural postoperative detachments
  - Enables suspension of the dura mater

##### > CLINICAL BENEFITS

- Reduced pain Aesthetics
- Improved quality of life

#### 7. CONTRAINDICATIONS

> List of contraindications likely to compromise the success of implantation :

- Acute or chronic, local or systemic infections
- Missing cover tissues
- Allergy to any of the ingredients
- Unstable fractures, destruction, or poor bone quality that may affect the stability of the devices

> Surgeons will need to assess each of their patients to determine the specific risks inherent in the procedure and the benefits for their patients.

#### 8. WARNINGS

> Prior to use, it is the responsibility of the surgeon to become acquainted with the implant procedures described in the documentation provided by EUROS. The surgeon is responsible for properly conducting the intervention and must ensure that the implanted device will not pose any risk of impingement with other implants already present in the patient's body.

> The manufacturer is not responsible for any complications arising from improper diagnosis, and/or incorrect operating techniques, limitations of treatment methods, or inadequate access.

> The surgeon is responsible for ensuring that the imaging data provided as input data for the product design is consistent with the patient's anatomy at the time of implanting the custom device.

> Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and proper implant placement, are important factors in the successful use of this device by the surgeon.

Patients with contraindications should be disqualified.

#### 9. INFORMATIONS FOR PATIENTS

> Patients and the parent(s) (or legal guardian(s)) should be informed that the safety and durability of the implant depends on the patient's behaviour.

> The patient should be informed of possible postoperative complications. The surgeon must inform the patient of the potential risks and adverse effects of artificial prostheses and have their consent regarding the proposed procedure.

> Upon discharge from the healthcare facility, the patient must have received instructions and warnings from the surgeon related to the treatment and therapy to be followed, as well as any restrictions on physical activities and sports. The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, but also of the need for periodic medical follow-up. In particular, it is recommended not to perform physical activities which would lead to deterioration of the cranial vault and to use protective items accordingly. The surgeon should inform the patient that it is not recommended to scratch the wound or scar after the procedure.

> The life of the implant may be less than the patient's life. The surgeon must warn the patient that the device cannot reproduce the flexibility, strength, reliability or durability of normal healthy bone, and that the device may need to be replaced in the future.

#### 10. PRECAUTIONS

> The fixation of the cranial vault is the surgeon's choice. Most often, the fixation is made by screws and the plate used for cranio-maxillofacial surgery. Suspension of the dura mater can be performed with surgical wire.

> Screwing in the fixation components and the fixation holes on the implant (if necessary) is to be carried out outside of the implant site. The length of any bone screws used must be less than the thickness of the custom-made cranial vault indicated on the label and on the operating support attached to the device. EUROS does not provide the anchoring system.

> Use caution when handling implants. Implants must not be scuffed or otherwise damaged. Damaged implants and implants with damaged packaging must not be used, and should be returned to EUROS.

> The implant should be stored in its original sealed packaging and in a clean, dry place. This packaging must not be exposed to direct sunlight, ionising radiation or extreme temperatures.

> This device is for single use only and should not be reused. Reuse of a single-use device that comes into contact with blood, bone, tissue, or other bodily fluids may result in injury to the patient or user. Risks associated with the reuse of a single-use device may include, but are not limited to, mechanical failure and transmission of infectious agents.

#### 11. ADVERSE EFFECTS

> The non-exhaustive list of adverse reactions includes:

- Infection
- Allergic reactions
- Venous thrombosis, pulmonary embolism
- Cardiovascular disorders
- Haematoma

> Treatment of certain adverse reactions may require additional surgery

#### 12. STERILITY AND PACKAGING

> The custom Crown cranial vault is delivered sterile (gamma irradiation) at a sterility assurance level (SAL) of 10<sup>-6</sup>.

> Check the integrity of each package before use.

> Do not use the device if the integrity of the sterile double bag is damaged (loss of vacuum, punctured bags, damaged sealing), or if the expiry date has been exceeded.

> Once the packaging is opened, the implant must never be resterilised. Any resterilisation may have adverse effects on the implant materials and cause premature failure.

#### 13. MRI INFORMATION

The safety and compatibility of the Crown device have not been evaluated in a magnetic resonance (MR) environment. Heating, migration, or image artefacts for the Crown device was not evaluated in a magnetic resonance (MR) environment. Nevertheless, the Crown implant is manufactured from non-ferromagnetic materials as mentioned in the "Description" section of these instructions for use.

#### 14. ELIMINATION

For any disposal following a storage, use, or ablation error, implants must follow the hospital waste disposal chain in accordance with environmental laws, strict hygiene rules, local regulations and procedures in force at the institution. In the event of a product returned to the manufacturer due to deficiency, the explant must be sent decontaminated.

#### 15. ADDITIONAL INFORMATION

EUROS is available to provide information on the Crown cranial vault to surgeons and medical personnel.

Reporting an incident:

Any healthcare professional who has a complaint or reason for dissatisfaction concerning the quality, identification, reliability, safety, efficacy and/or performance of the device, must inform EUROS. In the event of a serious incident or risk of a serious incident likely to result in or has resulted in the death or serious deterioration of a patient or user's health, EUROS must be notified immediately by telephone, fax or e-mail (qualite@euros.fr). Defective products must be returned to EUROS decontaminated.

For further information regarding the device, please contact EUROS or your distributor.



EURO SAS

Z.E. Athélie III - 13600 La Ciotat - France

Tel: +33 4 42 71 42 71 | Fax: +33 4 42 71 42 80

www.euros-orthopaedics.com

#### 16. MEANING OF SYMBOLS

	Sterilised by irradiation
	Do not reuse
	Do not resterilize
	Attention: Consult the instructions for use and the operating technique
	Do not use if the packaging is damaged
	Protect from moisture
	Keep away from sunlight
	Expiry date
	See the precautions for use
	Lot code
	Catalogue reference
	Manufacturer
	Fragile; handle carefully
	Manufacture date
	Temperature limit

#### 17. MEANING OF ABBREVIATIONS

Qty: quantity

Mat: material

fr: French

en: English

it: Italian

es: Spanish

tr: Turkish

de: German

zh: Chinese

The reference text is the French text.