# COSYX

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FR

NI.D2002A1 – 2023/07/07

EN

Instructions for use

IT

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ES

Nota de instrucciones

DE

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NL

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PT

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## 0. REVISION

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INSTRUCTIONS FOR USE  
FOR SURGEONS AND MEDICAL STAFF

## 0. REVISION

Version & date	Modifications	URL internet
A1 2023/07/07	Creation	eifu.euros.fr/NI.D200 2A.pdf

## 1. DESCRIPTION

- > EUROS hip prostheses are intended to reconstruct the coxofemoral joint.  
The COSYX acetabular reinforcement is intended for the reconstruction of the acetabulum and is used with cemented cup.
- > The materials used to manufacture these products satisfy current standard requirements (see the end of this leaflet).

## 2. INDICATIONS

The indications for use of total hip replacement prostheses include:

- Severely painful and/or disabled joint as a result of osteoarthritis, post-traumatic arthritis, rheumatoid arthritis or congenital hip dysplasia;
- Avascular necrosis of the femoral head;
- Acute traumatic fracture of the femoral head or neck;
- Failure of a previous hip surgery.

## 3. CONTRAINDICATIONS

Hip replacement is contraindicated in the following cases:

- Acute, systemic or chronic infection;
- Muscular, neurological or vascular deficiency of the affected limb;
- Pathologies that may compromise the functionality of the implant in any way;
- Bone destruction, or loss of bone characteristics that may compromise the stability of the implant;
- Allergy to any prosthesis component.

## 4. WARNINGS

- > The use of this product is reserved for surgeons with sufficient skills and experience in hip replacement surgery. The product must be used under the conditions specified by the manufacturer for the designated purpose.
- > Prior to use, surgeons shall be trained to the assembly

and installation techniques and shall have read EUROS documentation. Insertion must be carried out in accordance with the recognised surgical technique. The surgeon is responsible for ensuring that the operation proceeds smoothly and correctly.

- > Compliance with preoperative and intraoperative procedures, including knowledge of surgical techniques and the selection and proper positioning of devices, are important factors in the successful use of this device or system by the surgeon. Furthermore, the selection of a suitable patient together with the latter's cooperation are essential factors for a successful outcome.  
Patients who present contra-indications must be turned down.
- > The surgeon should take into consideration factors that may affect the performance of the device or system:
  - The patient's weight, activity level and occupation;
  - The patient's bone state (osteoporosis, tumors);
  - The patient's allergies;
  - The patient's diseases and disorders (infectious diseases, mental or neuromuscular disorders);
  - The patient's addiction (drug, alcoholism or nicotine).

The surgeon must assess the patient's physical and mental state and decide whether the latter is suited to having the device fitted, can stand up to the operation and follow a suitable rehabilitation program.

- > The manufacturer is not liable for any complications arising from incorrect diagnosis, incorrect choice of device, the incorrect assembly of components of the device and/or incorrect operating techniques, the limitations of treatment methods or inadequate aseptic conditions.
- > Reporting an incident:  
Any healthcare professional wishing to make a complaint or express their dissatisfaction in terms of the quality, identification, reliability, safety, efficacy and/or the performance of EUROS device must inform EUROS.  
In the case of a serious incident or risk of a serious incident liable to lead to or that may have led to death or to severe deterioration in the state of health of a patient or user, EUROS must be notified immediately by phone, fax or email ([complaints@euros.fr](mailto:complaints@euros.fr)). The competent authority of the Member State in which the user/or patient is established must be notified in case of a serious incident in relation to the device.  
For any additional information involving the devices, please contact EUROS or your distributor.  
EUROS SAS  
Z.E. Athélia III - 13600 La Ciotat - France

## 5. PRECAUTIONS

- > Caution should be exercised when handling the devices. The devices must not be scratched or damaged in any way whatsoever. Damaged devices and those with damaged packaging must not be used and should be returned to EUROS.
- > Devices must be kept in its original sealed packaging in a clean, dry location. Store the sterile products between 10 and 30°C. During transport, temperatures of up to 40°C for short periods can be tolerated. The packaging must not be directly exposed to the sun's rays or extreme temperatures.
- > Devices are for single use only and must not be reused. The reuse of a single use device that has come into contact with blood, bone, tissue or other body fluids may cause injury to the patient or user. The risks associated with the reuse of a single-use device may include, but are not limited to, mechanical failure and the transmission of infectious agents.
- > Devices must only be used with the instruments specified by EUROS in the reference surgical technique.
- > In the event of any disposal due to an error of storage or use, or to excision, devices should be disposed of following the disposal method for hospital waste in accordance with the procedures in force in the establishment.
- > Using all means available such as X-rays and templates, the surgeon must ensure that, by guaranteeing the concordance of scales between the X-rays and the templates, the implantation of the device will cause no risks of impingement with other implants already in the patient's body.

### **COSYX : Acetabular reinforcement for cemented cup**

- > Only to be used with bone defects. This Kerboul-type support ring must be adapted to the acetabulum. It is recommended to check that the posterior hook is correctly positioned in the obturator foramen and check the hold of the screws in the superior section.
- > The acetabulum fixation must be performed using  $\varnothing 4.5$  cancellous bone screws.

## 6. EXPECTED PERFORMANCE AND BENEFITS

- > Technical performances
  - Restoring the bearing capacity of the acetabulum:

- Restoration of a normal hip center ;
- Restoration of acetabular bone stock ;
- Enabling stable cemented fixation:
  - Screw fixation and hook shape;
- Being adapted to the acetabulum anatomy;
- Enabling reliable assembly with cemented acetabular cups:
  - Maintaining the fragmented graft to full integration ;
  - Resisting to the mechanical loads applied over its lifetime.
- Answering the implantation constraints:
  - Biocompatible ;
  - Being easily cleaned and sterilized ;
  - Being easily implantable and removable.
- > Clinical performances
  - Restoration of load bearing capacity of the acetabulum ;
  - Restoration of hip joint stability, kinematics and mobility.
- > Clinical benefits
  - Decreased pain ;
  - Increased quality-of-life ;
  - Restoring hip mobility.

## 7. ADVERSE EFFECTS / RESIDUAL RISKS

The list of potential adverse events associated with any hip arthroplasty includes:

- Excessive wear of the implant components secondary to impingement of components or damage of articular surfaces ;
- Instability ;
- Fracture, migration, loosening, dissociation, subsidence, protrusion, impingement, subluxation, or dislocation of the prosthesis or any of its components ;
- Possible detachment of the coating(s) on the femoral stem or acetabular shell components, potentially leading to increased debris particles ;
- Periprosthetic bone fracture, including femoral or acetabular perforation while seating the device ;
- Acetabular erosion or protrusion ;
- Osteolysis, geode ;
- Soft or surrounding tissues injuries ;
- Metal sensitivity reactions or other allergic/histological reactions to implant materials; inflammation, formation of erythema or oedema or implant rejection ;
- Superficial or deep infection, septic loosening ;

- Delayed wound healing ;
- Vascular damage resulting in blood loss and/or hematoma, potentially requiring transfusion ;
- Neurologic injury or neuropathy resulting in transient or permanent weakness, pain, and/or numbness ;
- Muscle or soft tissues injury ;
- Undesirable leg lengthening or shortening, claudication ;
- Periarticular calcification or heterotopic ossification, with or without impediment to joint mobility ;
- Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction ;
- Gastrointestinal or genitourinary complications ;
- Aggravation of other joint or back conditions due to positioning during surgery or neurological injury ;
- Decreased range of motion ;
- Squeaking (for ceramic-on-ceramic bearings) ;
- Intractable pain ;
- \*Intoxication, cellular damage, convulsion ;
- \*Cancer, tumor, immunological, and neurological disorders, genetic defects, non-heritable adverse effects in the progeny and/or an impairment of male or female reproductive functions or capacity.

Treatment of certain undesirable side effects could necessitate additional surgery.

*\*Note: these residual risks are related to the device biocompatibility. The biological risk assessment report demonstrated that the biological risk was controlled and considered acceptable for the EUROS products.*

## 8. INFORMATION TO BE PROVIDED TO THE PATIENT

- > The surgeon should inform the patient that the device's safety and durability will be affected by her/his weight, behaviour, namely the physical activities. Before leaving the care facility, the patient must receive instructions and warnings from the surgeon regarding the treatment and therapy to follow, and any restriction on physical and sporting activity or practices. The surgeon must inform the patient of all physical and psychological limitations inherent in the use of the device, but also of the need for regular medical follow-up. The lifespan of the implant may be shorter than the patient's lifetime or any particular given length of time. The surgeon must inform the patient that the device cannot and does not reproduce the flexibility, strength, reliability or

durability of normal healthy bone, and that the device may need to be replaced in the future.

- > MR environment: the surgeon must inform the patient to indicate that he is wearing an implant before an MRI examination.

An MRI safety evaluation based on a bibliographical analysis and a demonstration of equivalence has demonstrated the EUROS hip prostheses are MR-conditional. A patient with these devices can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0-Tesla ;
- Maximum spatial gradient field of 1500-Gauss/cm ;
- Normal operating mode only ;
- Whole-body-averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning for patient in a 3T MR system.

MR imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the EUROS hip prosthesis.

- > The surgeon must inform the patient to report any incident in relation to the devices to the surgeon.
- > The surgeon should inform the patient about the residual risks available in section 7 of this instruction for use.

## 9. CHOICE OF DEVICE / DEVICE COMPATIBILITY

- > The surgeon is responsible for selecting the devices according to the patient's morphology and the nature of the condition.

This choice may require preoperative planning. The instrumentation provided with the implant can also guide this choice thanks to a range of trial implants corresponding to the range of implants.

- > Information on device size is provided on the device packaging and on the device itself. The surgical team must check this information before the surgery.
- > In the case of devices designed to be implanted with cement or bone substitute, no incompatibility of use has been detected. The selection of the cement is the discretion of the surgeon.

Regarding devices designed to be implanted without cement, the manufacturer cannot, under any circumstances, be held responsible for implanting the device with cement.

## COMBINATION OF DEVICES

> EUROS devices must be combined only with EUROS devices. All combination with devices of other manufacturer is forbidden.

## 10. STERILITY

- > EUROS hip prostheses are supplied sterile (gamma irradiation) at a sterility assurance level (SAL) of  $10^{-6}$ .
- > Inspect every package before use. Do not use the component if the integrity of the double envelope ensuring the maintenance of sterility has been compromised (loss of vacuum, damaged envelopes, damaged seal area), or if the expiry date has elapsed.
- > Once the package has been opened, the device must never be re-sterilised. Any re-sterilisation may have adverse effects on the device materials and cause premature defects.
- > For instruments used to insert implants, see the relevant leaflet for information on sterilisation.

## 11. ADDITIONAL INFORMATION

The Summary of Safety and Clinical Performance is available upon request (email: [qualite@euros.fr](mailto:qualite@euros.fr) ; Tel : +33 4 42 71 42 71).

COSYX basic UDI-DI : 37004372H00197

## 12. MEANING OF THE SYMBOLS

List of symbols applicable for the implants covered by these instructions for use

	Sterilized using irradiation Double sterile barrier system
	Do not reuse
	Caution
	Batch code
	Keep dry
	Keep away from sunlight
	Storage temperature limitation
	Do not re-sterilise
	Use-by date
	Do not use if package is damaged and consult the instructions for use

	Catalogue number
	Manufacturer
	Date of manufacture
	Consult the instructions for use
	Medical device
	UDI

## 13. MEANING OF THE ABBREVIATIONS

Qty: quantity - fr: French - en: English - it: Italian - es: Spanish - de: German – nl : Dutch – pt : Portuguese  
SS: Stainless steel

## 14. MATERIAL

COSYX	100% X2CrNiMo 18-15-3 (ISO 5832-1 and ASTM F139)
-------	--

## MATERIAL CHEMICAL COMPOSITION

Stainless steel X2CrNiMo 18-15-3 (ISO 5832-1 and ASTM F139)

%C	%Si	%Mn	%P	%S	%N
≤0.030	≤0.75	≤2.0	≤0.025	≤0.010	≤0.10
%Cr	%Mo	%Ni	%Cu	%Co	%Fe
17.0 to 19.0	2.25 to 3.0	13.0 to 15.0	≤0.50	≤0.10	Balance

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